


BULBWORKS

SPECIALTY LIGHTBULBS AND SOCKETS

PRICE SHEET & CROSS REFERENCE GUIDE

July, 2015



Replacement Lamps For: Ophthalmic Equipment & Precision Optical Instruments

Replacement Lamps For American Optical;
Bausch & Lomb; Burton; Copeland; Edgemaster;
Haag Streit; Humphrey; Marco; Mentor; Neitz; Nikon;
Reichert; Thorpe; Topcon; Welch Allyn; Zeiss & Many Others

BULBWORKS, INC.
P.O. BOX 586 · SUCCASUNNA, NJ 07876

CALL: 1-800-334-2852 · 1-973-584-7171
FAX: 1-800-335-2858 · 1-973-584-0300

E-mail: bulbwork@bulbworks.com · Visit Us On The Web At: www.bulbworks.com



Replacement Lamps For:

Ophthalmic Equipment & Precision Optical Instruments

Price Sheet & Cross Reference Guide

July, 2015

(Good Through July 31, 2016)

Welcome...

Price Sheet.....	1
Illustrations (Figure Table).....	5
Cross Reference Guide.....	7
Terms & Conditions.....	11

Pay Less, Expect More...

Just because our lamps cost less doesn't mean you'll get less! See for yourself why Bulbworks is fast catching the attention of Dealers, Equipment Manufacturers and Specialty Lamp Buyers all over the country. We will quickly become a valuable asset to your business.

Bulbworks offers only the highest quality specialty lamps for use in Ophthalmic Equipment, Microscopes, Inspection Equipment, Comparators, Optical, Medical, & Scientific Equipment. We feature lamps from Ushio, Philips, Osram/Sylvania, GE, Welch Allyn, Wiko, and many other top manufacturers –

The same brands used by the equipment manufacturers, but you'll save big!

You can buy with confidence - All lamps are fully warranted!
We will immediately replace any defective lamp!

We Go Out Of Our Way For You...

At Bulbworks, customer service is not a department - *it's our philosophy!* Expect knowledgeable, courteous, and memorable service. For reliability, we can't be beat! With Toll Free Phone and Fax Numbers, we're convenient, economical, and we ship *fast!*

BULBWORKS, INC. • P.O. BOX 586 • SUCCASUNNA, NJ 07876

CALL TOLL FREE: 800-334-2852 • (973-584-7171)

FAX TOLL FREE: 800-335-2858 • (973-584-0300)

E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com

Replacement Lamps For Ophthalmic Equipment & Precision Optical Instruments

Prices Good Through July 31, 2016

Save On The Bulbs You Need! Call For Any Item Not Listed!

<u>Industry Stand.Code</u>	<u>Description</u>	<u>Bulbworks P/N</u>	<u>Price Each</u>	<u>Fig. #</u>
<i>ANSI Coded Lamp Types</i>				
BRL	12V,50W G6.35 Base	BW.BRL	\$8.05	3
BRL-OSRAM 64610	12V,50W Osram Brand	BW.BRL-OSR	8.65	3
BZW	120V,50W SC Bayonet	BW.BZW	15.95	54
CAX-130V (Import)	130V,50W DC Bayonet	BW.CAX-130V	7.75	52
CAX-130V (GE Brand)	130V,50W DC Bayonet	BW.CAX-130V-GE	9.90	52
CDD (Import)	120V,100W SC Bayonet	BW.CDD	3.90	54
CDD-GE (GE Brand)	120V,100W SC Bayonet	BW.CDD-GE	CALL	54
DNE	120V,150W MR16 Vented	BW.DNE	33.12	70
EFR	15V,150W MR16	BW.EFR	23.15	69
ESB 7388 (Philips)	6V,20W G4 Bi-pin	BW.ESB	6.51	2
ESB-OSRAM 64250	6V,20W G4 Bi-pin	BW.ESB-OSR	6.51	2

Miniature Lamp Types

50	7.5V, .22A Miniature Screw	BW.50	\$1.62	42
238	6.3V, .50A Miniature Bayonet	BW.238	3.65	51
327	28V, .04A T1 3/4 Midget Flange	BW.327	2.65	23
328	6V, .20A SC Midget Flange	BW.328	2.65	23
328-GE (GE Brand)	6V, .20A SC Midget Flange	BW.328-GE	3.20	23
605	6.15V, .50A Min.Screw	BW.605	1.62	58
787	6V,10W G4 Bi-pin Base	BW.787	13.58	1
1096	6V,4.5A DC Pref.Base	BW.1096	17.65	9
1460 (Import)	6.5V,2.75A DC Pref. Base	BW.1460	11.15	10
1460-GE (GE Brand)	6.5V,2.75A DC Pref. Base	BW.1460-GE	CALL	10
1468 (Import)	6V,4.5A DC Pref. Base	BW.1468	11.40	11
1468-GE (GE Brand)	6V,4.5A DC Pref. Base	BW.1468-GE	14.95	11
1489	6.5V,2.75A SC Bayonet	BW.1489	12.90	47
1493 (Import)	6.5V,2.75A DC Bayonet	BW.1493	7.45	44
1493-GE (GE Brand)	6.5V,2.75A DC Bayonet	BW.1493-GE	8.62	44
1507	12.5V,3.0A Pref. Base	BW.1507	11.80	-
1680X	6V,4.10A SC Bayonet	BW.1680X	10.95	45
6251	6V,5W SC Bayonet	BW.6251	5.15	42

Wattage Coded Lamp Types

6S6 DC-120V	120V,6W DC Bayonet	BW.6S6.DC-120V	\$2.54	48
10S6-120V	120V,10W Cand.Screw	BW.10S6-120V	3.35	59
10S11N	120V,10W Interm.Screw Clear	BW.10S11N	4.50	64
10S11N/IF	130V,10W Interm.Screw Frost	BW.10S11N.CW130	1.00	66
15S11/IF	120V,15W Cand.Base Frost	BW.15S11.IF	14.61	62

Wattage Coded Lamp Types Are Continued On the Next Page

Bulbworks, Inc. CALL: 1-800-334-2852 (973-584-7171) • FAX: 1-800-335-2858 (973-584-0300)
E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com

Replacement Lamps For Ophthalmic Equipment & Precision Optical Instruments
Prices Good Through July 31, 2016

Save On The Bulbs You Need! Call For Any Item Not Listed!

<u>Industry Stand.Code</u>	<u>Description</u>	<u>Bulbworks P/N</u>	<u>Price Each</u>	<u>Fig. #</u>
<i>Wattage Coded Lamp Types Are Continued From The Previous Page</i>				
15S11/13 (Import)	120V,15W Cand.Base Clear	BW.15S11.13	\$5.95	62
15S11/13-GE (GE Brand)	120V,15W Cand.Base Clear	BW.15S11.13-GE	7.80	62
15S11/14	115-120V,15W Cand.Base	BW.15S11.14	9.99	62
15S11/87CL	120V,15W Intermediate Clear	BW.15S11.87 CL	7.62	64
15S11/87FR	120V,15W Intermediate Frosted	BW.15S11.87FR	8.15	65
15S11/102CF	120V,15W Med. Screw Frost	BW.15S11.102CF	5.50	68
15S11/102CL-GE	120V,15W Med. Screw Clear	BW.15S11102CLGE	5.37	67
15S11/102CL	120V,15W Med. Screw Clear	BW.15S11.102CL	4.86	67
15T7/C	120V,15W Cand. Base	BW.15T7C	2.77	60
15T7.DC-120V	120V,15W DC Bayonet	BW.15T7.DC-120V	3.39	50
25S11/N	115-125V,25W Interm.Base	BW.25S11N	7.75	65
25S11/2C	120V,25W Cand.Base Clear	BW.25S11.2C	9.10	62
25S11/5C	120V,25W Cand.Base Frost	BW.25S11.5C	9.10	63
25T8/SCP	6V,25W SC Pref.Base	BW.25T8.SCP	37.10	8

Equip.Mfrs.Part Numbers (Number Sequences Are Grouped by Manufacturer)

1003	12.8V,.94A SC Bayonet	BW.1003	\$3.47	46
1008	6V,4.5A G12½ 9 Hole Base	BW.1008	25.05	15
2007-6	6V,4.5A S11 Spec.Pref.Base	BW.510-250	CALL	-
2011-B (New Perimeter)	8V,12W DC Bayonet	BW.2011-B	CALL	47
2011-S (Old Perimeter)	6V,12W Screw Base	BW.2011-S	CALL	-
3007	120V,15W Intermediate Base	BW.15T7N	6.75	61
4007	115-125V,25W Interm.Base	BW.25S11N	7.75	65
4010	12V,4.2A Special Base	BW.45070	CALL	72
76907	6V,15-18W DC Bayonet Base	BW.76907	12.80	43
77904	6V,30W DC Bayonet Base	BW.77904	14.36	49
77910	6V,15W Pref.Flange Base	BW.77910	14.90	16
78079	6V,15W DC Triang.Pref.Base	BW.78079	17.05	7
88052	3V,0.8A Spec.Threaded Base	BW.88052	CALL	-
10080	6V,2.4W Frosted E10 Base	BW.10080	11.37	58
11496-802	2.7V,0.45A Spec.Screw Base	BW.11496	32.92	18
11530	6V,15W SC Bayonet Base	BW.11530	16.35	45
12230	6V,5W DC Bayonet Base	BW.12230	11.55	40
14800	6V,33W SC Spec.Triang.Base	BW.14800	29.10	12
18200	6V,33W Pref.Ring Base	BW.18200	CALL	-
20030	12V,30W G6.35 Bi-Pin Base	BW.20030	39.75	4
20570	6V,3W E10 Screw Base	BW.20570	9.95	57
24400	6V,15W Pref.Flange Base	BW.24400	38.65	16

Equipment Manufacturer Part Numbers Are Continued On The Next Page

Bulbworks, Inc. CALL: 1-800-334-2852 (973-584-7171) • FAX: 1-800-335-2858 (973-584-0300)

E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com

Replacement Lamps For Ophthalmic Equipment & Precision Optical Instruments
Prices Good Through July 31, 2016

Save On The Bulbs You Need! Call For Any Item Not Listed!

<u>Industry Stand.Code</u>	<u>Description</u>	<u>Bulbworks P/N</u>	<u>Price Each</u>	<u>Fig. #</u>
<i>Equip.Mfrs.' Part Numbers Are Continued From The Previous Page (Number Sequences Are Grouped by Manufacturer)</i>				
25500	6V,33W DC Spec.Flange Base	BW.40300-25500	CALL	8
40340-20700	6V,20W Prefocused Base	BW.44740-20400	37.60	71
40350-42110	6V,20MA G4 Bi-Pin	BW.7628	3.47	-
40770	6V,.2A Spec.Single Pin Base	BW.40770	24.75	25
41366	6V,33W DC Spec.Flange Base	BW.40300-25500	CALL	-
42409-30060	7V,35W DC Bayonet Base	BW.42409-30060	19.85	48
80020	6V,10W Miniature Flange	BW.80020	CALL	-
L-04	4V,3.6W Special Screw Base	BW.L-04	36.50	-
L-09	12.5V,60W SC Pref.Base	BW.L-09	46.40	14
L-13	4V,0.3W T1 Spec.Screw Base	BW.L-13	CALL	-
71-71-37	7V,35W S8 DC Pref.Base	BW.71-71-37	28.15	11
71-71-55	2.5V,0.3A Spec.Threaded Base	BW.71-71-55	42.50	27
71-75-62	2.5V,0.38A Spec.Thread Base	BW.71-75-62	11.95	35
510-250	6V,4.5A S11 Spec.Pref. Base	BW.510-250	CALL	-
HS366	6V,4.5A S11 9 Hole Base	BW.HS366	31.50	15
900-930	6V,4.5A S11 9 Hole Base	BW.HS366	31.50	15
OQ-53Z	6V,25W Triangle Flange Base	BW.OQ-53Z	29.80	6
OQ-58Z	6V,30W Pref.Flange Base	BW.OQ-58Z	34.20	13
M34	6V,20W G4 Bi-pin	BW.M34	14.70	2
45070	12V,50W Spec. 9 Hole Screw Base	BW.45070	CALL	72
450012	12V,20W G4 T3 Bi-Pin Base	BW.450012	CALL	2
JC12V-50W	12V,50W GY6.35 Base	BW.JC12V-50W	11.34	4
80A	60V,7W Ring Fluorescent	BW.80A	141.25	-
64260	12V,30W Special 2 Pin Base	BW.64260	35.40	5
64425	12V,20W G4 Bi-pin Base	BW.64425	7.65	3

Welch Allyn Lamp Types

WA00200	2.5V Vacuum	BW.WA00200	\$24.85	29
	2.5V Vacuum	BW.00200 (Equiv.Brand)	18.35	29
WA00600	2.5V Vacuum	BW.WA00600	14.85	35
WA00900	2.5V Vacuum	BW.WA00900	26.95	37
WA01100	2.5V Vacuum	BW.WA01100	CALL	22
WA01800	2.5V Vacuum	BW.WA01800	13.20	39
WA02200	2.5V Vacuum	BW.WA02200	13.20	21
WA02500	6V Vacuum	BW.WA02500	24.42	56
WA02600	6V Halogen	BW.WA02600	46.40	24
WA02900	2.5V Vacuum	BW.WA02900	CALL	-
WA03000	3.5V Halogen	BW.WA03000	47.87	32
	3.5V Halogen	BW.03000 (Equiv.Brand)	25.65	32

Welch Allyn Lamp Types Are Continued On The Next Page

Bulbworks, Inc. CALL: 1-800-334-2852 (973-584-7171) • FAX: 1-800-335-2858 (973-584-0300)

E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com

Replacement Lamps For Ophthalmic Equipment & Precision Optical Instruments
Prices Good Through July 31, 2016

Save On The Bulbs You Need! Call For Any Item Not Listed!

<u>Industry Stand.Code</u>	<u>Description</u>	<u>Bulbworks P/N</u>	<u>Price Each</u>	<u>Fig. #</u>
<i>Welch Allyn Lamp Types Are Continued From The Previous Page</i>				
WA03100	3.5V Halogen	BW.WA03100	\$51.30	30
	3.5V Halogen	BW.03100 (Equiv.Brand)	25.65	30
WA03300	2.5V Halogen	BW.WA03300	CALL	-
WA03400	2.5V Halogen	BW.WA03400	45.73	31
	2.5V Halogen	BW.03400 (Equiv.Brand)	25.65	31
WA03600	2.5V Vacuum	BW.WA03600	12.35	55
WA03700	3.5V Halogen	BW.WA03700	41.40	38
	3.5V Halogen	BW.03700 (Equiv.Brand)	26.95	38
WA03900	2.5V Halogen	BW.WA03900	57.95	34
	2.5V Halogen	BW.03900 (Equiv.Brand)	14.05	34
WA04100	14.5V Halogen	BW.WA04100	53.05	26
	14.5V Halogen	BW.04100 (Equiv.Brand)	14.25	26
WA04400	2.5V Halogen	BW.WA04400	15.00	33
	2.5V Halogen	BW.04400 (Equiv.Brand)	10.50	33
WA04500	3.5V Halogen	BW.WA04500	44.75	32
	2.5V Halogen	BW.04600 (Equiv.Brand)	CALL	27
WA04700	2.5V Vacuum	BW.WA04700	27.50	28
	2.5V Vacuum	BW.04700 (Equiv.Brand)	11.90	28
WA04800	2.5V Vacuum	BW.WA04800	12.40	28
	2.5V Vacuum	BW.04800 (Equiv.Brand)	11.90	28
WA04900	3.5V Halogen	BW.WA04900	64.95	33
	3.5V Halogen	BW.04900 (Equiv.Brand)	21.45	33
WA06000	3.5V Halogen	BW.WA06000	44.95	29
	3.5V Halogen	BW.06000 (Equiv.Brand)	21.45	29
WA06100	14.5V Halogen	BW.WA06100	15.65	27
	14.5V Halogen	BW.06100 (Equiv.Brand)	21.45	27
WA06200	3.5V Halogen	BW.WA06200	CALL	20
WA07200	2.5V Halogen	BW.WA07200	41.70	34
WA07600	2.5V Halogen	BW.WA07600	CALL	-
WA07800	6V Vacuum	BW.WA07800	46.63	19
	6V Vacuum	BW.07800 (Equiv.Brand)	21.50	19
WA08100	12V,20W T3 G4 Bi-pin	BW.64425	7.65	3
WA08200	3.5V Halogen	BW.WA08200	29.65	36
	3.5V Halogen	BW.08200 (Equiv.Brand)	25.65	36
WA08300	3.5V Halogen	BW.WA08300	CALL	36
	3.5V Halogen	BW.08300 (Equiv.Brand)	CALL	36
WA09500	24W Solarc Lamp	BW.WA09500	475.00	73
WA72000	2.5V Rechargeable Battery	BW.WA72000	CALL	74
WA72200	3.5V Rechargeable Battery	BW.WA72200	CALL	75
WA72600	2.5V Blue Rechargeable Battery	BW.WA72600	CALL	-

Long Life Halogen Replacement For CAX-130V

RN13050	130V,50W DC Bayonet	BW.RN13050-USH	\$35.10	53
---------	---------------------	----------------	---------	----

Many Other Lamp Types Are Available. Please Call For Any Item Not Listed!
Bulbworks, Inc. CALL: 1-800-334-2852 (973-584-7171) • FAX: 1-800-335-2858 (973-584-0300)
E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com

Lamp Illustrations May Represent A Variety Of Different Lamps




Fig. 1 Fig. 2 Fig. 3 Fig. 4 Fig. 5 Fig. 6 Fig. 7 Fig. 8



Fig. 9 Fig. 10 Fig. 11 Fig. 12 Fig. 13 Fig. 14



Fig. 15 Fig. 16 Fig. 17 Fig. 18 Fig. 19 Fig. 20 Fig. 21



Fig. 22 Fig. 23 Fig. 24 Fig. 25 Fig. 26 Fig. 27 Fig. 28 Fig. 29 Fig. 30



Fig. 31 Fig. 32 Fig. 33 Fig. 34 Fig. 35 Fig. 36 Fig. 37 Fig. 38 Fig. 39

Lamp Illustrations May Represent A Variety Of Different Lamps



Fig. 40



Fig. 41



Fig. 42



Fig. 43



Fig. 44



Fig. 45



Fig. 46



Fig. 47



Fig. 48



Fig. 49



Fig. 50



Fig. 51



Fig. 52



Fig. 53



Fig. 54



Fig. 55



Fig. 56



Fig. 57



Fig. 58



Fig. 59



Fig. 60



Fig. 61



Fig. 62



Fig. 63



Fig. 64



Fig. 65



Fig. 66



Fig. 67



Fig. 68



Fig. 69



Fig. 70



Fig. 71




Fig. 72



Fig. 73



Fig. 74

Fig. 75

Catalog Effective Date Of July, 2015
 No Part Of This Catalog May Be Reproduced Without Permission.
 Pictures Are For Reference Only And Are Not To Scale.
BULBWORKS, INC. · P.O. BOX 586 · SUCCASUNNA, NJ 07876
 Call: 1-800-334- (BULB) 2852 · 1-973-584-7171 Fax: 1-800-335-2858 · 1-973-584-0300
 E-mail: bulbwork@bulbworks.com Visit Us On The Web At: www.bulbworks.com



Replacement Lamps For Ophthalmic Equipment & Precision Optical Instruments
Cross Reference Guide (Prices Good Through July 31, 2016)

Here's Just A Partial List...

<u>Manufacturer</u>	<u>Mfr.Part No.</u>	<u>Description</u>	<u>Bulbworks P/N</u>	<u>Price</u> <u>Ea</u>	<u>Fig.</u> <u>No.</u>
<u>AIT</u>					
630 Layout Blocker		6V,4.10A SC Bayonet	BW.1680X	\$10.95	45
Essilor Posicentron		6V,5W SC Bayonet	BW.6251	5.15	42

<u>American Optical</u>					
	10407	120V,6W DC Bayonet	BW.6S6.DC-120V	\$2.54	48
Projector	10409	120V,10W Cand.Screw	BW.10S6-120V	3.35	59
Projector	10454	6.15V,.50A Min.Screw	BW.605	1.62	58
Projector	10462	6.3V,.50A Min.Bayonet	BW.238	3.65	51
	11063	6.5V,2.75A DC Pref.Base	BW.1460	11.15	10
	11143	6V,20W G4 Bi-Pin Halogen	BW.ESB	6.51	2
Lensometer	11163	130V,10W Frosted Int.Screw	BW.10S11N.CW130	1.00	66
Projectors 1215;1217	11176	130V,50W Clear	BW.CAX-130V	7.75	52
		130V,50W Clear (GE Brand)	BW.CAX-130V-GE	9.90	52
Lensometer 11210	11211	120V,15W Frosted Int. Screw	BW.15S11.87FR	8.15	65
	11360	120V,15W Frosted Int. Screw	BW.15S11.87FR	8.15	65
	11581	6V,4.5A DC Prefocused Base	BW.1468	11.40	11
Fixation	11583	7.5V,.22A Miniature Screw	BW.50	1.62	42
Proj-O-Chart 11800	11807	6V,4.5A DC Prefocused Base	BW.1096	17.65	9
Tonotarget	12419	6V,.20A SC Midget Flange	BW.328	2.65	23
Lensometer	12603	120V,15W Clear Int.Screw	BW.15S11.87 CL	7.62	64
Lensometer	12604	120V,10W Int. Screw	BW.10S11N	4.50	64
Proj-O-Marker		6.5V,2.75A DC Bayonet	BW.1493	7.45	44
	12610	120V,15W DC Bayonet	BW.15T7.DC-120V	3.39	50

<u>Bausch & Lomb</u>					
Vertometer 70 Frosted	21-26-62	120V,25W Frosted Cand. Screw	BW.25S11.5C	\$9.10	63
Lensometer 70 Frosted	21-65-18	120V,25W Frosted Cand. Screw	BW.25S11.5C	9.10	63
Lensometer 70 Clear	21-65-21	120V,25W Clear Cand.Screw	BW.25S11.2C	9.10	62
Lensometer		120V,15W Clear Cand.Screw	BW.15S11.13	5.95	62
Lensometer		120V,15W Frosted Cand. Base	BW.15S11.IF	8.80	62
Projector	71-34-40	115-125V,50W SC Bayonet	BW.BZW	15.95	54
Slit Lamp OSL2000;PSL3000	71-61-86	12V,50W 9 Hole Spec.Screw Base	BW.45070	CALL	72
Projector Acuity Chart	71-71-37	7V,35W S8 DC Pref. Base	BW.71-71-37	28.15	11
Perimeter Auto Plot	71-71-41	6.5V,2.75A SC Bayonet	BW.1489	12.90	47
Slit Lamp	71-71-50	6V,25W SC Pref. Bayonet	BW.25T8.SCP	37.10	8
Retinoscope	71-71-55	2.5V,0.3A Spec.Threaded Base	BW.71-71-55	42.50	27
Projector	71-71-77	115-125V,50W SC Bayonet	BW.BZW	15.95	54
Projector Compact Acuity	71-71-78	120V,100W T8 SC Bayonet	BW.CDD	3.90	54
		120V,100W T8 SC Bayonet	BW.CDD-GE	CALL	54
Perimeter	71-71-80	12.5V,3.0A Prefocused Base	BW.1507	11.80	-
Keratometer 71-21-35	71-71-84	115-120V,15W Cand.Screw	BW.15S11.14	9.99	62
Retinoscope	71-75-62	2.5V,0.38A Spec. Threaded Base	BW.71-75-62	19.95	35

Replacement Lamps For Ophthalmic Equipment & Precision Optical Instruments

Cross Reference Guide (Prices Good Through July 31, 2016)

Here's Just A Partial List...

<u>Manufacturer</u>	<u>Mfr.Part No.</u>	<u>Description</u>	<u>Bulbworks P/N</u>	<u>Price Ea</u>	<u>Fig. No.</u>
<u>Burton</u>					
20/20 Lensometer		120V,15W Med. Scrw.Frosted	BW.15S11.102CF	\$5.50	68
Slit Lamp 800;850;1000		6V,4.5A G12½ 9 Hole Base	BW.1008	25.05	15
Slit Lamp 2000		12V,50W GY6.35 Base	BW.JC12V-50W	11.34	4
Keratometer 2040		60V, 7W Ring Fluorescent	BW.80A	141.25	-
Projector 5000		130V,50W DC Bayonet	BW.CAX-130V	7.75	52
		130V,50W DC Bayonet	BW.CAX-130V-GE	9.90	52

<u>Copeland</u>					
Retinoscope	71-71-55	2.5V,0.3W Spec.Threaded Base	BW.71-71-55	\$42.50	27

<u>Edgemaster</u>					
Blocker		6.5V,2.75A DC Bayonet	BW.1493	\$7.45	44

<u>Haag-Streit</u>					
Perimeter	510-250	6V,4.5A Spec.Prefocused Base	BW.510-250	CALL	-
Slit Lamp 900-930	366	6V,4.5A S11 9 Hole Base	BW.HS366	\$31.50	15

<u>Humphrey</u>					
Visual Field Perimeter		12.8V,.94A SC Bayonet	BW.1003	\$3.47	46
Visual Field Analyzer		12V,50W G6.35 Bi-Pin Base	BW.BRL	8.05	3

<u>Marco</u>					
Models EL2B; 3B; 5B; 5G; 2BB-R; 2B-Photo; 5R; 5 Photo; Slit Lamp					
	1008	6V,4.5A G12½ 9 Hole Base	BW.1008	\$25.05	15
	1036	120V,150W MR16 Vented 2 Pln	BW.DNE	33.12	70
Slit Lamp G2; G4; G20	1045	12V,30W Special G4 Base	BW.64260	35.40	5
Perimeter	2007-6	6V,4.5A S11 Special Pref.Base	BW.510-250	CALL	-
Perimeter White Chart (Newer)	2011-B	8V,12W DC Bayonet	BW.2011-B	CALL	47
Perimeter White Chart (Older)	2011-S	6V,12W Screw Base	BW.2011-S	CALL	-
	2020	15V,150W MR16	BW.EFR	23.15	69
Keratometer	3006	120V,15W Clear Cand.Base	BW.15T7C	2.77	60
	3007	120V,15W Intermed.Screw Base	BW.15T7N	6.75	61
Keratometer	3008	120V,15W Clear Cand.Base	BW.15T7C	2.77	60
101 Lensometer	4007	115-125V,25W Intermd.Scrw.Frost.	BW.25S11N	7.75	65
Perimeter	4008	6V,4.5A DC Prefocused Base	BW.1468	11.40	11
Radioscope 5005	5010	6.5V,2.75A DC Bayonet	BW.1493	7.45	44
Projector	6014	130V,50W DC Bayonet	BW.CAX-130V	7.75	52
		130V,50W DC Bayonet	BW.CAX-130V-GE	9.90	52

<u>Mentor</u>					
Slit Lamp	1008	6V,4.5A G12½ 9 Hole Base	BW.1008	\$25.05	15
Slit Lamp Models S10; G10; SH10; GH12	22-4010	12V,50W 9 Hole Spec.Screw Base	BW.45070	CALL	72

Replacement Lamps For Ophthalmic Equipment & Precision Optical Instruments

Cross Reference Guide (Prices Good Through July 31, 2016)

Here's Just A Partial List...

<u>Manufacturer</u>	<u>Mfr.Part No.</u>	<u>Description</u>	<u>Bulbworks P/N</u>	<u>Price Ea</u>	<u>Fig. No.</u>
<u>Neitz</u>					
Retinoscope R1 Streak	L-04	4V,3.6W Special Screw Base	BW.L-04	\$36.50	-
Slit Lamp SL-J (Main)	L-09	12.5V,60W SC Pref.Base	BW.L-09	46.40	14
Slit Lamp SLJ; SLCX (Fix.)	L-13	4V,,3W T1 Special Screw Base	BW.L-13	CALL	-
Lensometer LM-P2	L-19	120V,15W Med.Screw Frosted	BW.15S11.102CF	5.50	68
<u>Nikon</u>					
Lensometer PL1	77904	6V,30W DC Bayonet	BW.77904	\$14.36	49
Lensometer Black Zoom PL2	78079	6V,15W DC Prefocused Base	BW.78079	17.05	7
Slit Lamp Models FS2; FS3;CS2; NS1	79761	12V,30W 2 Pin Prefocused Base	BW.64260	35.40	5
Retinoscope	88052	3V,0.8A Special Threaded Base	BW.88052	CALL	-
<u>Reichert</u>					
Ophthalmoscope	11063	6.5V,2.75A DC Pref. Base	BW.1460	\$11.15	10
		6.5V,2.75A DC Pref. Base	BW.1460-GE	CALL	10
Project-O-Chart	11082	130V,50W DC Bayonet	BW.CAX-130V	7.75	52
		130V,50W DC Bayonet	BW.CAX-130V-GE	9.90	52
Projector 11083; 11075	11143	6V,20W G4 Bi-Pin Halogen	BW.ESB	6.51	2
Fulvue Retinoscope	11496-802	2.7V,0.45A Spec.Screw Base	BW.11496	32.92	18
Ophthalmoscope	11581	6V,4.5A DC Prefocused	BW.1468	11.40	11
<u>Thorpe</u>					
Slit Lamp		6V,25W SC Prefocused Base	BW.25T8.SCP	\$37.10	8
<u>Topcon</u>					
Ophthalmoscope OM3 & 4	10080	6V,2.4A Frosted E10 Screw Base	BW.10080	\$11.37	58
Lensometer LMP5 & 6	10550	6V,15-18W DC Bayonet Base	BW.76907	12.80	43
Ophthalmoscope OM1	11530	6V,15W SC Bayonet Base	BW.11530	16.35	45
Perimeter SBP-20	12230	6V,5W DC Bayonet Base	BW.12230	11.55	40
Perimeter SBP-11 Main	14800	6V,33W SC Spec.Triangular Base	BW.14800	29.10	12
Perimeter SBP-20 Main	18200	6V,33W Pref. Ring Base	BW.18200	CALL	-
Projector ACP5/6S/6R/7R	20030	12V,30W G6.35 Bi-Pin Base	BW.20030	17.40	4
Perimeter SBP-11	20570	6V,3W E10 Screw Base	BW.20570	9.95	57
Slit Lamp Models 2D; 1E; 2E; 4E; 7E	20700	6V,20W Prefocused Base	BW.44740-20400	37.60	16
Slit Lamp SL-2	24400	6V,15W Prefocused Base	BW.24400	38.65	-
Slit Lamp SL1; 1A; 1D	25500	6V,33W DC Special Flange	BW.40300-25500	CALL	26
Slit Lamp Models SL3D; SL5D; SL3E; SL6E (Main Bulb)	25600	6V,4.5A G12½ 9 Hole Base	BW.1008	25.05	15
Main Bulb	25600	6V,4.5A G12½ 9 Hole Base	BW.1008	25.05	15
Fixation Bulb	40310-40770	6V,20MA Single Pin Base	BW.40770	24.75	25
Fixation Bulb	40350-42110	6V,20MA G4 Bi-Pin Base	BW.7628	3.47	-
Indirect Ophthalmoscope	41500-30130	6V,10W G4 Bi-Pin Base	BW.787	13.58	1
Lensometer LMP3 & 4	30060	7V,35W DC Bayonet Base	BW.42409-30060	19.85	-

Lamps for Topcon Are Continued On the Next Page

Replacement Lamps For Ophthalmic Equipment & Precision Optical Instruments

Cross Reference Guide (Prices Good Through July 31, 2016)

Here's Just A Partial List...

<u>Manufacturer</u>	<u>Mfr.Part No.</u>	<u>Description</u>	<u>Bulbworks P/N</u>	<u>Price</u> <u>Ea</u>	<u>Fig.</u> <u>No.</u>
Topcon Lamps Are Continued From Previous Page					
<u>Topcon</u>					
Projector Models ACP4D; ACP-1; CP-3; CP-4D					
	30060	7V,35W DC Bayonet	BW.42409-30060	\$19.85	-
Projector CP5D	31600	130V,50W DC Bayonet	BW.CAX-130V	7.75	52
		130V,50W DC Bayonet	BW.CAX-130V-GE	9.90	52
Lensometer	T-23 Clear	120V,15W Clear, Med. Screw	BW.15S11.102CL	4.86	67
Lensometer	T-23 Frost	120V,15W Frosted Med.Screw	BW.15S11.102CF	5.50	68
Lensometer	51010	120V,15W Frosted Med.Screw	BW.15S11.102CF	5.50	68
Sceenscope	51010	120V,15W Frosted Med.Screw	BW.15S11.102CF	5.50	68
Keratometer OM2	51010	120V,15W Frosted Med.Screw	BW.15S11.102CF	5.50	68
Lensometer LM6E	51010	120V,15W Frosted Med.Screw	BW.15S11.102CF	5.50	68
Ophthalmoscope ID-10	80020	6V,10A Miniature Flange	BW.80020	CALL	-

<u>Welch Allyn</u>					
Otoscope	WA00200	2.5V Vacuum	BW.WA00200*	\$24.85	29
Ophthalmoscope/Retinoscope	WA00900	2.5V Vacuum	BW.WA00900	26.95	37
Retinoscope	WA01800	2.5V Vacuum	BW.WA01800	13.20	39
Headlight	WA02500	6V Vacuum	BW.WA02500	24.42	56
Headlight	WA02600	6V Halogen	BW.WA02600	46.40	24
Standard Ophthalmoscope	WA03000	3.5V Halogen	BW.WA03000*	47.87	32
Transilluminator	WA03100	3.5V Halogen	BW.WA03100*	51.30	30
Ophthalmoscope	WA03300	2.5V Halogen	BW.WA03300	CALL	-
Otoscope	WA03400	2.5V Halogen	BW.WA03400*	45.73	31
18100 Retinoscope	WA03700	3.5V Halogen	BW.WA03700*	41.40	38
Ophthalmoscope	WA03900	2.5V Halogen	BW.WA03900*	57.95	34
Ophthalmoscope	WA04400	2.5V Halogen	BW.WA04400	15.00	33
Retinoscope	WA04500	3.5V Halogen	BW.WA04500	44.75	32
Ophthalmoscope	WA04900	3.5V Halogen	BW.WA04900	64.95	33
Streak Retinoscope	WA07200	2.5V Halogen	BW.WA07200	CALL	34
Retinoscope	WA08200	3.5V Halogen	BW.WA08200*	29.65	36
Retinoscope	WA08300	3.5V Halogen	BW.WA08300*	CALL	36
Rechargeable Battery	WA72200	3.5V Rechargeable Battery	BW.WA72200	CALL	75

*Also Available in Generic Brand - See Price Section

<u>Woodlyn</u>					
Slit Lamp HR1;HR3;HR5;STR	45070	12V,50W 9 Hole Spec.Screw Base	BW.45070	CALL	72
Slit Lamp	45500	6V,20W Spec.G4 Base	BW.44740-20400	37.60	71
Keratometer		115-120V,15W Cand.Screw Base	BW.15S11.14	9.99	62
Slit Lamp	45570	6V,20W Spec.G4 Base	BW.44740-20400	37.60	71
Slit Lamp	450012	12V,20W T3 G4 Bi-Pin Base	BW.450012	CALL	2

<u>Zeiss</u>					
Slit Lamp 100/16;125/16	39-01-53	6V,25W Spec.Triang.Flange	BW.OQ-53Z	\$29.80	6
Slit Lamp 30SL/M	3801-20-7040	12V,30W Spec.G4 Base	BW.64260	35.40	5

Many Other Lamp Types Are Available. Please Call For Any Item Not Listed!

Bulbworks, Inc. CALL: 1-800-334-2852 (973-584-7171) • FAX: 1-800-335-2858 (973-584-0300)

E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com

Getting To Know Us... (Terms & Conditions)

ORDERING IS EASY...

Call Us TOLL FREE • 800-334-BULB(2852) • Our business hours are 8:30 AM to 5:00 PM Eastern Time Monday through Thursday and 8:30 AM to 4:00 PM on Friday. Outside of business hours? Use our after hours message machine and your call will be returned during the next business day or you may...

Fax Us ON OUR 24 HOUR TOLL FREE FAX LINE • 800-335-2858 • All orders received by fax are acknowledged promptly either by phone or return fax.

Order Online • www.bulbworks.com • All orders placed online will be acknowledged promptly by return e-mail.

No Minimum Order • We appreciate your business no matter how small or large your order may be. Orders with a merchandise value of under \$10.00 will be subject to a Small Order Handling Charge of \$3.00 (In addition to Merchandise and Shipping).

Fast Shipment • Orders received by 3:00 PM Eastern Time will be **shipped the same day** that we receive them. Orders are normally shipped by UPS but all express services are available. Bulbworks carefully inspects all lamps prior to shipment and we use the best packaging methods available.

Terms • We accept MasterCard, Visa, Discover, & American Express. Our terms for *approved* open accounts are 2% 10 Days, Net 30, freight prepaid and added to the invoice. We would be glad to open an account for your convenience. Just ask!



Order With Confidence • All items carry a full warranty. Bulbworks sells only top quality brand name products - the same brands used by the equipment manufacturers. Should product failure, breakage, or defect occur, contact us for immediate replacement!

Original Manufacturer's Packaging • Each lamp is shipped in its original packaging from the lamp manufacturer – Eiko, Philips, GE, Ushio, Osram - Sylvania, Welch Allyn, and others. Part numbers are clearly labeled and are easy to read to eliminate confusion at reorder time. Helpful equipment reference information appears on every carton.

All Prices Are Subject To Change Without Notice. Some Prices May Be For Stock On Hand Only.
No Part Of This Catalog May Be Reproduced Without Permission.



P.O. BOX 586 SUCCASUNNA, NJ 07876
CALL: 800-334-BULB (2852) • (973-584-7171)
FAX: 800-335-2858 • (973-584-0300)

E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com



**Need Something Besides Ophthalmic Bulbs? Stop Looking & Call Us!
We Have What You Need! We Stock All Types Of
Bulbs, Sockets, and Batteries Including:**

- **MICROSCOPE & SCIENTIFIC BULBS** Replacement Lamps For American Optical, Bausch & Lomb, Jena, Leitz, Nikon, Olympus, Swift, Unitron, Zeiss, And Others. Lamps For Biological, Biomedical, Educational, Forensic, Industrial, Metallurgical, Research, And Other Applications. Fluorescent And LED Microscope Illuminators.
- **MEDICAL BULBS** Replacement Lamps For Operatory, Surgical, Clinical, Fiber Optic, Laser, Analytical, Veterinary, X-Ray Viewing, And Other Applications. Lamps For All Types Of Medical Equipment.
- **INSPECTION** Replacement Lamps For Optical Comparators Including Baty Gagemaster, J & L, Micro Vu, Mitutoyo, Nikon, OGP, Scherr Tumico, Vermont Precision, And Others. Replacement Lamps for Microscopes And Other Inspection Devices.
- **GRAPHIC ARTS BULBS** Replacement Lamps For AB Dick, Agfa, Amergraph, Brown, Burgess, Douthitt, Itek, Kenro, NuArc, Olec, Sandmar, Screen Theimer, 3M, And More.
- **OPTICAL AND OPHTHALMIC BULBS** Replacement Lamps For American Optical, Bausch & Lomb, Copeland, Haag Streit, Marco, Mentor, Nikon, Reichert, Topcon, Welch Allyn, And Others. Lamps For Optical And Ophthalmic Equipment Including Lensometers, Ophthalmoscopes, Projectors, Retinoscopes, Slit Lamps, Tonometers, And Others.
- **GENERAL AND SPECIALTY LIGHTING** Replacement Lamps And LEDs For Color Enhancement And Correction, Commercial, Decorative, Energy Saving, General Service, Long Life, Office, Problem Solving, Rough Service, And Other Applications.
- **GERMICIDAL, BLACK LIGHT, AND OTHER UV BULBS** Replacement Lamps For Air And Water Purification. For Use In Hospitals, Laboratories, Manufacturing, Food Handling, And Many Other Applications.
- **PROJECTOR AND AUDIO VISUAL BULBS** Replacement Lamps For Overhead Projectors, Slide And Movie Projectors, Microfilm and Microfiche Readers, Micrographic And Other Applications.
- **SOCKETS AND LAMP HOLDERS** Replacement Sockets For All Lamp Types Including Halogen, Incandescent, Vacuum, & Others. Sockets From Bender & Wirth, Buhl, Gilway, Christiana Industries, Cole Hersee, Gilway, Sylvania, Ushio, And Others.
- **BATTERIES** Batteries For All Types Of Devices. Alkaline, Lithium, Nickel Metal Hydride, Heavy Duty, Zinc Air, Rechargeable, Industrial, Hearing Aid, Photo, Electronic, Medical, Keyless Entry, Watch and Others. All Sizes Available.

We Supply All Specialty Lamp Types Including: Fluorescent, Incandescent, Halogen, Mercury Vapor, Metal Halide, Miniature, Quartz, Xenon, And Other Specialty Bulbs.

BULBWORKS Supplies Only Top Quality Brand Names - The Same Brands Used By The Equipment Manufacturers. Lamps From Eiko, GE, Hosobuchi, Narva, Nikon, Opti Quip, Osram Sylvania, Philips, Ushio, Welch Allyn, ORC, Voltarc, Zeiss, And Many Other Specialty Lamp Manufacturers

Call BULBWORKS for Your Next Replacement Lamp Requirement!

BULBWORKS, INC. • P.O. BOX 586 • SUCCASUNNA, NJ 07876

CALL: 800-334-BULB(2852) • (973-584-7171)

FAX: 800-335-2858 • (973-584-0300)

E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com