



LED And Fluorescent Microscope Illuminators



BW.HX-1
(Fluorescent)

110V, 8W, Daylight, 6400 Degrees K, Ring Shape, Fluorescent Ring Light Assembly For Microscopes, Self Contained Unit Includes Bulb, Power Supply And Connecting Cable, Maximum Connecting Caliber Is 60mm, Includes 1 Additional Bulb
Uses Replacement Bulb BW.92ADL

Dealer Net Price: \$85.95



BW.HX-2
(Fluorescent)

The Fluorescent Version Of This Item BW.HX-2 Has Been Discontinued. The Nearest Replacement Is The LED Version BW.HXD-7. Replacement Fluorescent Bulbs For HX-2 Remain Available As Our Part Number BW.FCB92D

Discontinued
See BW.HXD-7



BW.HX-4
(Fluorescent)

120V, 8W, Daylight, 6400 Degrees K, Ring Shape, Fluorescent Ring Light For Microscopes, Self Contained Unit Includes Bulb, Power Supply, Switch And Connecting Cable, Maximum Connecting Caliber Is 55mm
Uses Replacement Bulb BW.FCB80A

Dealer Net Price: \$79.50



BW.HXD-1
(LED)

120V, 6W, Daylight, 6000 Degrees K, Ring Shape, LED Ring Light For Microscopes, 60 LEDs, 96-265 VAC, Dimmable, Includes Plastic Power Supply With Controller, Maximum Connecting Caliber Is 60mm, Lamp Head Is ABS Plastic And Aluminum

Dealer Net Price: \$124.95



BW.HXD-2
(LED)

120V, 12W, 3200 Degrees K, Gooseneck Microscope Illuminator With Self Contained Power Supply, 2 Tubes, LED, 6W Per Tube, 300mm (11.8") Gooseneck Length, 92-220 VAC

Dealer Net Price: \$131.95



BW.HXD-2B
(LED)

120V, 16W, Daylight, 6000 Degrees K, Gooseneck Microscope Illuminator With Self Contained Power Supply And Heavy Duty Metal Base, 500mm (19.7") Gooseneck Length, Dual Head, LED, Each Head 8 Watts, Dimmable, 50-60 HZ, 85-265 VAC

Dealer Net Price: \$168.95



BW.HXD-3
(LED)

120V, 6W, Daylight, 6000 Degrees K, Ring Shape, Ring Light For Microscope, Includes Ring Light, Heavy Duty Metal Power Supply With Controller And Connector, 60 LEDs (6 Watt), Dimmable, 96-220 VAC, Maximum Connecting Caliber Is 60mm, Lamp Head Is Sturdy Plastic And Metal

Dealer Net Price: \$148.95



BW.HXD-4
(LED)

120V, 3W, Daylight, 6500-7000 Degrees K, Ring Light For Microscope, Self Contained Unit Includes Power Supply with Line Cord, Switch and Ring Light, Dimmable, Maximum Connecting Caliber Is 60mm, Durable Plastic Housing

Dealer Net Price: \$69.95



BW.HXD-5
(LED)

120V, 6W, Daylight, 6000 Degrees K, Ring Shape, Ring Light For Microscope, 60 LEDs (6 Watt), 96-220 VAC, Dimmable, Maximum Connecting Caliber Is 60mm, Plastic Power Supply Case, Lamp Head Is Durable Plastic And Aluminum

Dealer Net Price: \$124.95



BW.HXD-6
(LED)

12W, Daylight, 6000 Degrees K, Ring Shape, Heavy Duty Ring Light For Microscope, Multi Level Lighting With Up To 12 Watts (80 LEDs), 5 Levels Of Lighting Intensity - Full or Additional Sub Area Lighting 1/2, 1/4, 1/8, Max. Connecting Caliber Is 60mm, Heavy Duty Power Supply & Lamp Head

Dealer Net Price: \$225.00



BW.HXD-7
(LED)

120V, 6W, Daylight, 6000 Degrees K, Ring Shape, Ring Light For Microscope, Self Contained Unit Includes Power Supply, Switch and Ring Light, Dimmable, Maximum Connecting Caliber Is 60mm, Durable Plastic Housing

Dealer Net Price: \$124.95



BW.HXD-9
(LED)

120V, 16W, Gooseneck Microscope Illuminator For Mounting Directly On Scope, Including Mounting Brackets, Dual LED Light Heads Each 3-8 Watts, Dimmable, 590mm (23.2") Gooseneck Length, 50-60 HZ, 85-265 VAC

Dealer Net Price: Call



BW.HXP-295
(LED)

120V, 4W, Circular Shape, LED Stage Illuminator For Microscope, 95mm Diameter x 4.88mm Thick, 4 Watt Assembly Includes Power Supply And Connector, Light Is For Bottom Lighting, Provides Transmitted Light, Dimmable, 90-265 VAC

Dealer Net Price: \$122.95



BW.92ADL
(Fluorescent)

8 Watt Replacement Lamp for HX-1 Microscope Illuminator

Dealer Net Price: \$64.95



BW.FCB92D
(Fluorescent)

10 Watt Replacement Lamp for HX-2 Microscope Illuminator

Dealer Net Price: \$115.40



BW.FCB80A
(Fluorescent)

8 Watt Replacement Lamp for HX-4 Microscope Illuminator

Dealer Net Price: \$42.75