


BULBWORKS

SPECIALTY LIGHTBULBS AND SOCKETS

Replacement Lamps For DENTAL EQUIPMENT

CONSUMER

July, 2016



Price Sheet & Cross Reference Guide

Replacement Lamps For Operatory, Curing, Fiberoptic,
Handpiece, Examination, Video, X-ray Viewing & More

Bulbworks, Inc.

P.O. Box 586, Succasunna, NJ 07876

Call: 800-334-2852 · (973-584-7171)

Fax: 800-335-2858 · (973-584-0300)

E-mail: bulbwork@bulbworks.com

Visit Us Online: www.bulbworks.com



Replacement Lamps For Dental Equipment Price Sheet & Cross Reference Guide

Replacement Lamps For Operatory, Curing, Fiberoptic,
Handpiece, Examination, Video, X-Ray Viewing & More

Welcome...

Price Sheet.....	Page 2
Full Spectrum/Color Corrected Lamps.....	Page 5
Pictures (Figure Table).....	Page 6
Equipment To Lamp Cross Reference Guide.....	Page 8
Distributor/Industry Number Cross Reference Guide.....	Page 15
Dental Waste Recycling Products.....	Page 17
Terms & Conditions.....	Page 18

Pay Less, Expect More...

Just because our lamps cost less doesn't mean you'll get less! Learn why Bulbworks has earned the business of Dental Dealers and Equipment Manufacturers all over the country. We will quickly become a valuable asset to *your* business.

Bulbworks offers only the highest quality specialty lamps for use in Dental Equipment. We feature lamps from Ushio, Philips, Osram/Sylvania, GE, Eiko, and many other top manufacturers - the same brands used by the equipment manufacturers, but you'll save big!

You can buy with confidence - All lamps carry a full warranty!
We will immediately replace any defective lamp!

We Go Out Of Our Way For You...

At Bulbworks, customer service is not a department - *it's our philosophy!* Expect knowledgeable, courteous, and memorable service. For reliability, we can't be beat! With Toll Free Phone and Fax Numbers, we're convenient, economical, and *we ship fast!*

Bulbworks, Inc. • P.O. Box 586 • Succasunna, NJ 07876
CALL: 800-334-2852 • (973-584-7171) • TOLL-FREE 24-HOUR FAX: 800-335-2858 (973-584-0300)
E-mail: bulbwork@bulbworks.com • Visit Us Online: www.bulbworks.com

**Replacement Lamps For Dental Equipment
Dealer Net Price Sheet, July, 2016 (Good Through July 31, 2017)
Operatory, Curing, Fiber Optic, Handpiece, Video**

Bulb	Description	Bulbworks Part Number	Efos No.	Retail Price*	Fig.
AL1824	24W Solarc Lamp	BW.WA09500 Welch Allyn Brand	---	\$515.00	8
AL1824	24W Solarc Lamp	BW.09500-USH Ushio Brand	---	409.50	8
AL1824	24W Solarc Lamp	BW.09500 Import Brand	---	205.50	8
BHB	120V,250W	BW.BHB	---	41.95	7
BU7012-GEN1	MK DENT QC5016 Coupler	BW.BU7012-GEN1 1 Pack Generic	---	38.60	42
BU7012-GEN3	MK DENT QC5016 Coupler	BW.BU7012-GEN3 3 Pack Generic	---	106.50	42
DNE	120V,150W	BW.DNE	3019	33.12	6
DNE USHIO	120V,150W	BW.DNE-USH Ushio Brand	3019	57.85	6
DNF	21V,150W	BW.DNF	3023	28.65	6
DNF USHIO	21V,150W	BW.DNF-USH Ushio Brand	3023	31.51	6
DP1020	WA Chair Site Mirror Bulb	BW.DP1020	---	48.95**	36
DP1520	Denlite Portable Mirror Bulb	BW.DP1520	---	36.45**	36
DZE	24V,150W	BW.DZE Osram Brand	---	19.10	17
EFM	8V,50W	BW.EFM	3040	21.85	2
EFN	12V,75W	BW.EFN	3035	23.15	1
EFN OSRAM	12V,75W	BW.EFN-OSR Osram Brand	3035	24.90	1
EFP	12V,100W	BW.EFP	3020	23.15	1
EFP OSRAM	12V,100W	BW.EFP-OSR Osram Brand	3020	25.45	2
EFR	15V,150W	BW.EFR	3014	23.15	1
EFR OSRAM	15V,150W	BW.EFR-OSR Osram Brand	3014	25.60	1
EHJ	24V,250W	BW.EHJ	3038	12.75	14
EHT	120V,250W	BW.EHT	4073	20.15	33
EJA	21V,150W	BW.EJA	3018	18.15	1
EJA OSRAM	21V,150W	BW.EJA-OSR Osram Brand	3018	23.15	1
EJA PHILIPS	21V,150W	BW.EJA-PHIL Philips Brand	3018	24.95	1
EJL	24V,200W	BW.EJL	3017	21.85	1
EJL OSRAM	24V,200W	BW.EJL-OSR Osram Brand	3017	24.05	1
EJL USHIO	24V,200W	BW.EJL-USH Ushio Brand	3017	24.95	1
EJM	21V,150W	BW.EJM	3011	26.65	1
EJM OSRAM	21V,150W	BW.EJM-OSR Osram Brand	3011	27.50	1
EJV	21V,150W	BW.EJV	3022	23.65	1
EJV EIKO	21V,150W	BW.EJV-EIKO Eiko Brand	3022	24.59	1
EJV OSRAM	21V,150W	BW.EJV-OSR Osram Brand	3022	25.45	1
EJV USHIO	21V,150W	BW.EJV-USH Ushio Brand	3022	26.37	1
EKE	21V,150W	BW.EKE	3027	16.26	1
EKE GE	21V,150W	BW.EKE-GE GE Brand	3027	20.79	1
EKE USHIO	21V,150W	BW.EKE-USH Ushio Brand	3027	18.60	1
ELC	24V,250W	BW.ELC	---	18.23	2
ELZ	21V,150W	BW.ELZ	3024	32.65	6
ENL	12V,50W	BW.ENL	3036	19.37	1
ENZ	30V,50W	BW.ENZ	4084	28.65	1
EPN	12V,35W	BW.EPN	---	38.15	1
EPV	14.5V,90W	BW.EPV	4080	24.85	1
EPX	14.5V,90W	BW.EPX	---	31.65	1
ESD	120V,150W	BW.ESD	3028	31.80	1
ESD USHIO	120V,150W	BW.ESD-USH Ushio Brand	3028	43.35	1
ESJ	82V,85W	BW.ESJ	---	36.52	1
ESP	120V,150W	BW.ESP	---	23.40	19
ESX	12V,20W	BW.ESX	4085	10.25	1
ESY	100V,150W	BW.ESY	3034	28.47	13
ETJ	120V,250W	BW.ETJ	3015	24.97	1
ETJ OSRAM	120V,250W	BW.ETJ-OSR Osram Brand	3015	26.45	1
EVA	12V,100W	BW.EVA	3059	14.15	12
EVA OSRAM	12V,100W	BW.EVA-OSR Osram Brand	3059	16.50	12
EZK	120V,150W	BW.JCR120V150WB	3008	31.75	1
FCS	24V,150W	BW.FCS	4070	7.69	14
FCS OSRAM	24V,150W	BW.FCS-OSR Osram Brand	4070	9.12	14
FDS/DZE	24V,150W	BW.FDS.DZE	4074	17.05	17
FDS/DZE	24V,150W	BW.FDS.DZE-OSR	4074	18.89	17
FDS/DZE	24V,150W	BW.FDS.DZE-USH	4074	19.65	17
FDT	12V,100W	BW.FDT	4081	14.73	18
FTD	12V,20W	BW.FTD	---	9.18	5

**Price Is For Stock On Hand Only

For Equipment To Lamp Information, Please Consult The Bulbworks, Inc. Cross Reference Guide

*Price shown here is the Bulbworks Suggested Retail Price. The Equipment Manufacturer's List Price for the same item may vary widely from what is indicated here.

Many Other Lamp Types Are Available! Please Call For Any Lamp Not Listed!

Bulbworks, Inc. CALL: 800-334-2852 (973-584-7171) • FAX: 800-335-2858 (973-584-0300)

E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com

Replacement Lamps For Dental Equipment
Dealer Net Price Sheet, July, 2016 (Good Through July 31, 2017)
Operatory, Curing, Fiber Optic, Handpiece, Video

Bulb	Description	Bulbworks Part Number	Efos No.	Retail Price*	Fig.
JA12V-55WD H3	12V,55W	BW.JA12V-55WD	4071	\$7.95	28
JA12V-100WD	12V,100W	BW.JA12V-100WD	---	11.40	28
JA24V-60WDBT	24V,60W	BW.JA24V60WDBT	---	29.50	29
JA24V-70WD	24V,70W	BW.JA24V-70WD	---	14.90	28
JC12V-50W	12V,50W	BW.JC12V-50W	3032	11.34	12
JC12V-50WG1	12V,50W	BW.JC12V-50WG1	---	13.05	12
JCD24V-60WBA	24V,60W	BW.JCD24V-60WBA	---	31.89	19
JCD24V-75W	24V,75W	BW.JCD24V-75W	---	32.25	19
JCD120V-150WB	120V,150W	BW.JCD120V150WB	---	26.42	12
JCR120V-150WB	120V,150W	BW.JCR120V150WB	3008	31.75	1
JCRM12V-50W	12V,50W	BW.JCRM12V-50W	---	33.19	3
JCRM12V-75W	12V,75W	BW.JCRM12V-75W	3016	30.15	4
JCRM12V-100W	12V,100W	BW.JCRM12V-100W	3021	35.60	4
M10P003	14W/9.5W MR11 Solarc Lamp	BW.M10P003	---	133.70	7
M21E001	60V,21W MR11 Welch Allyn Brand	BW.M21E001	---	372.00	9
M21E001-IMP	60V,21W MR11 Import Brand	BW.M21E001-IMP	---	190.50	9
M34	6V,20W	BW.M34	---	14.90	11
MHR150N	95V,150W	BW.MHR150N	---	315.00	10
PELTON & CRANE LAMP TOOL	To Install Lamp	BW.PC TOOL	4078	16.15	31
Q75CL/DC-28V	28V,75W	BW.Q75CL.DC-28V	4079	20.10	19
Q150T4/CL25V	25V,150W	BW.Q150T4.CL25V Ushio Brand	4076	19.63	30
Q150T4/CL25V IMP	25V,150W	BW.Q150T4CL25IM Import Brand	4076	13.52	30
SABLE 1	Access L & Rotomax L Handpiece	BW.SABLE-1	---	Call	36
SABLE 5	Access L & Rotomax L Handpiece	BW.SABLE-5	---	Call	36
THE71	15V,150W	BW.EFR	---	23.15	2
VRB-GEN1	Vector VKS6 Coupler Bulb	BW.VRB-GEN1 1 Pack Generic	---	Call	41
VRB-GEN3	Vector VKS6 Coupler Bulb	BW.VRB-GEN3 3 Pack Generic	---	Call	41
WA09500	24W Solarc Lamp	BW.WA09500 Welch Allyn Brand	---	515.00	8
WA09500	24W Solarc Lamp	BW.09500-USH Ushio Brand	---	409.50	8
WA09500	24W Solarc Lamp	BW.09500 Import Brand	---	205.50	8
WA09800	21W Metal Halide	BW.WA09800	---	Call	8
Y900529	NSK Handpiece Bulb	BW.Y900529-GEN1 Generic Brand -1 Pack	---	Call	41
Y900529	NSK Handpiece Bulb	BW.Y900529-GEN3 Generic Brand - 3 Pack	---	Call	41
3S6-120V	120V,3W	BW.3S6-120V	---	2.70	27
25A-120V	120V,25W	BW.25A-120V	---	3.27	21
125T10/P	120V,125W	BW.125T10.P	4072	32.25	20
300SB/IF	120V,300W	BW.300SB.IF	---	30.72	26
1222 LED	Kinetic Handpiece - LED Version - Rated 3.5V	BW.1222-LED Can Be Used To Replace Standard 1222 Bulb	---	68.25	44
13165	14V,35W	BW.13165 Philips Brand	3026	41.20	3
13165 IMPORT	14V,35W	BW.13165-IMP Import Brand	3026	20.75	3
13298	10V,52W	BW.3367-5267**	3009	52.11	3
13865	12V,75W	BW.13865	---	42.55	4
14552	12V,75W	BW.14552	---	49.55	3
21237	12V,80W	BW.21237	---	36.80	4
21237 EIKO	12V,80W	BW.21237-EIKO Eiko Brand	---	54.15	4
22076 LED	W & H Handpiece - LED Version - Rated 3.5V	BW.22076-LED Can Be Used To Replace Standard 22076 Bulb	---	66.95	45
60986	2 Pack Star Quick Connect Swivel	BW.60986	---	90.65	32
64115	12V,10W	BW.64115	---	15.10	55
64156	24V,70W	BW.64156	---	15.77	28
64440	12V,50W	BW.64440	---	8.30	12
64450AX	12V,75W	BW.64450AX	---	10.50	15
64451	24V,60W No Black Tip	BW.64451	---	32.80	28

**Same As Caulk/Dentsply #644751 - also fits the 3M XL1000 operated at different voltage.

For Equipment To Lamp Information, Please Consult The Bulbworks, Inc. Cross Reference Guide

*Price shown here is the Bulbworks Suggested Retail Price. The Equipment Manufacturer's List Price for the same item may vary widely from what is indicated here.

Bulbworks, Inc. CALL: 800-334-2852 (973-584-7171) • FAX: 800-335-2858 (973-584-0300)

E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com

**Replacement Lamps For Dental Equipment
Dealer Net Price Sheet, July, 2016 (Good Through July 31, 2017)
Operatory, Curing, Fiber Optic, Handpiece, Video**

Bulb	Description	Bulbworks Part Number	Efos No.	Retail Price*	Fig.
64602	12V,50W	BW.64602	---	\$26.45	12
64613	12V,75W	BW.64613-IMP	3010	24.85	4
64617	12V,75W	BW.64617	3016	32.23	4
64617S	12V,75W	BW.64617S	---	34.35	4
64624	12V,100W	BW.64624	3021	34.91	4
64638	24V,100W	BW.64638	---	24.20	14
64638	24V,100W	BW.64638-IMP	---	16.52	14
95775	41V,200W	BW.95775	---	40.85	1
561398	13V,100W	BW.561398	---	77.95	3
650064	2 Pack Midwest Power Optic For In-Hose Swivel	BW.650064-GEN2 Generic Brand – 2 Pack	---	61.40	38
650197	3 Pack Midwest Power Optic Handpiece Bulb	BW.650197 Midwest Brand	3025	Call	35
790257	Midwest Stylus Handpiece	BW.790257GEN1	---	63.95	36
790257 LED	Midwest Stylus Handpiece – LED Version	BW.790257-LED Can Be Used To Replace Standard 790257 Bulb	---	66.95	46
553-3881	2 Pack Kavo Lux Coupler Bulb – 3.3V,4W	BW.553-3881 Kavo Brand – 2 Pack	---	166.75	42
553-3881	1 Pack Kavo Lux Coupler Bulb – 3.3V,4W	BW.5533881-GEN1 Generic Brand – 1 Pack	---	38.60	42
553-3881	3 Pack Kavo Lux Coupler Bulb – 3.3V,4W	BW.5533881-GEN3 Generic Brand – 3 Pack	---	106.50	42
553-3881 LED	Kavo Lux Coupler Bulb – LED Version - Rated 3.3V, 4W	BW.553-3881-LED Can Be Used To Replace Standard 553-3881 Bulb	---	66.95	47
0.03523010	ADEC RA24 Roto-Quick Coupler	BW.0.03523010	---	86.75	43
041-117-00	ADEC Excellence (Older Style)	BW.041-117-00	---	58.20	39
041-317-00	ADEC "Newer" Style 3.5V	BW.041-317-00	3031	53.25	39
041-317-00 LED	ADEC "Newer" Style – LED Version – Rated 3.5V	BW.041-31700LED Can Be Used To Replace Standard 041-317-00 Bulb	---	66.95	48
041-709-00	17V,95W	BW.041-709-00	---	25.95	16
041-709-00	17V,95W	BW.041-709-00IM	---	19.75	16
3367-5267	10V,52W	BW.3367-5267**	3009	52.11	3
42001-02121	Dabi Atlante	BW.42001021GEN1 Generic Brand – 1 Pack	---	35.85	42
42001-02121	Dabi Atlante	BW.42001021GEN3 Generic Brand – 3 Pack	---	96.75	42
70201102004	13V,75W	BW.70201102004	---	177.50	2

Fluorescents For X-Ray Viewing; Magnifying Lights, Curing, Work Lights

Bulb	Description	Bulbworks Part Number	Efos No.	Retail Price*	Fig.
CF9DS.E.841	9W Comp. Fluor. 4100K 4 Pin Base	BW.CF9DS.E.841	---	\$17.95	60
F6T5/CW	6W 9" Cool White 5/8" Diam.	BW.F6T5.CW	3050	7.45	61
F6T5/D	6W 9" Daylight 5/8" Diam.	BW.F6T5.D	3048	8.25	61
F8T5/CW	8W 12" Cool White 5/8" Diam.	BW.F8T5.CW	3051	7.85	61
F8T5/D	8W 12" Daylight 5/8" Diam.	BW.F8T5.D	3049	8.60	61
F10T4/CW	10W 12 1/8" Cool White	BW.F10T4.CW	---	29.70	61
F13T5/CW	13W 21" Cool White 5/8" Diam.	BW.F13T5.CW	3052	10.30	61
F13T5/D	13W 21" Daylight 5/8" Diam.	BW.F13T5.D	---	11.10	61
F14T8/CW	14W 15" CW 1" Diam.	BW.F14T8.CW	3055	8.80	62
F14T8/D	14W 15" Daylight 1" Diam.	BW.F14T8.D	3056	10.30	62
F14T12/CW	14W 15"CW 1 1/2" Diam.	BW.F14T12.CW	3053	14.98	63
F14T12/D	14W 15" Daylight 1 1/2" Diam.	BW.F14T12.D	3054	15.35	63
F15T8/CW	15W 18" Cool White 1" Diam.	BW.F15T8.CW	3057	10.30	62
F15T8/D	15W 18" Daylight 1" Diam.	BW.F15T8.D	3058	10.60	62
FC8T9/CW	22W 8" Diameter Ring	BW.FC8T9.CW	8023	15.95	58

For Equipment To Lamp Information, Please Consult The Bulbworks, Inc. Cross Reference Guide

*Price shown here is the Bulbworks Suggested Retail Price. The Equipment Manufacturer's List Price for the same item may vary widely from what is indicated here.

**Same As Caulk/Dentsply #644751 - also fits the 3M XL1000 operated at different voltage.

**Bulbworks, Inc. CALL: 800-334-2852 (973-584-7171) • FAX: 800-335-2858 (973-584-0300)
E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com**

**Replacement Lamps For Dental Equipment
Dealer Net Price Sheet, July, 2016 (Good Through July 31, 2017)
X-Ray Duplicator Lamp**

Bulb	Description	Bulbworks Part Number	Efos No.	Retail Price*	Fig.
F8T5/BLB	8W 12" Blacklight Blue 5/8" Diam.	BW.F8T5.BLB	3049B	\$18.25	64

X-Ray Darkroom Lamps; Safelights

Bulb	Description	Bulbworks Part Number	Efos No.	Retail Price*	Fig.
7C7-ORANGE	120V,7W Orange, Candelabra	BW.7C7-ORANGE	---	\$5.26	25
7C7-RED	120V,7W Red, Candelabra	BW.7C7-RED	---	5.26	25
7½S11/CR	7½W Red, Med. Screw	BW.7.5S11.CR	---	6.15	24
15A/R	120-130V,15W Red Medium Screw Base	BW.15A.R	---	14.70	23
15A15-120V	120V,15W	BW.15A15-120V	---	5.10	22
25A/R	120V,25W	BW.25A.R	---	9.59	23
F15T8/R	15W 18" Red Fluor. 1" Diam.	BW.F15T8.R	---	31.89	65

Indicator, Panel, Communicator, And Miniature Lamps

Bulb	Description	Bulbworks Part Number	Efos No.	Retail Price*	Fig.
40	6.3V, .15A T3 1/4 Miniature Screw	BW.40	---	\$1.62	-
47	6.3V, .15A Miniature Bayonet	BW.47	---	1.62	54
85	28V, .04A Wedge Base	BW.85	---	2.95	40
88	6.8V, 1.91A S8 DC Bayonet	BW.88	---	3.10	57
88 GE	6.8V, 1.91A S8 DC Bayonet	BW.88-GE GE Brand	---	6.10	57
120PSB5	120V, .025A Slide Base	BW.120PSB5	---	3.77	49
336	14V, .08A Midget Groove	BW.336	---	1.62	52
376	28V, .06A T1 3/4 Midget Flange	BW.376	---	2.82	53
381	6.3V, .20A T1 3/4 Midget Flange	BW.381	---	3.40	53
387	28V, .04A Midget Flange	BW.387	---	2.03	53
387 GE	28V, .04A Midget Flange	BW.387-GE GE Brand	---	3.10	53
787	6V, 10W G4 Base	BW.787	---	14.02	34
789	12V, 14W T2 3/4 G4 Base	BW.789	---	18.35	34
1073	12.8V, 1.80A SC Bayonet	BW.1073	---	1.62	56
1819	28V, .04A T3 1/4 Miniature Bay	BW.1819	---	3.65	
7382	14V, .08A T1 3/4 Bi-Pin	BW.7382	---	3.07	37
462420	24V, 20MA Slide Lamp	BW.462420	---	3.50	50
552450	24V, 50MA Slide Lamp	BW.552450	---	4.95	51
HPR40	6V, .67A T3 SC Mini Flange	BW.HPR40	---	2.75	55
HPR53	4V, 3.4W T3 SC Mini Flange	BW.HPR53	---	5.78	55

Full Spectrum Lighting/Color Corrected

Full Spectrum Color Corrected Fluorescent Lamps • For The General Lighting Of Dental Offices And Laboratories

Benefits Of Full Spectrum Color Corrected Lamps Include:

- Accurate Shade And Color Matching! • Diagnostic Precision! • Visual Accuracy! • Greater Efficiency Of Personnel! • Increased Eye Comfort! • Office Area Is Made More Pleasant And Invigorating! • Helps People To Work And Feel Better!

Bulb	Description	Bulbworks Part Number	Retail Price*	Fig.
F15T8/DSGN50	15W 18" Color Corrected 1" Diam.	BW.F15T8.DSGN50	\$19.50	62
F32T8/850	32W 48" Length 1" Diameter Color Corrected	BW.F32T8.850 25 Pieces Per Case. Brand Can Vary	6.95	62
FO32T8/850	32W 48" Length, 1" Diameter Color Corrected	BW.FO32.850 30 Pieces Per Case, Sylvania or SLI Brand	8.50	62
F40T12.C50 Replaces Standard F40T12CW	48" Length, 1½" Diameter	BW.F40T12.C50 30 Pieces Per Case Only	11.95	63
FB40T12EXCU6 (Excella) Replaces Standard FB40T12CW	6" Leg Spacing U Shaped	BW.FB40T12EXCU6 Limited Availability	Call	-

Call For Any Size Or Type Of Fluorescent Lamp Not Listed. Energy Saving Versions Are Also Available.

Color Corrected Fluorescents Fit Most Standard Fluorescent Light Fixtures With No Modification Necessary.

Because of the fragile nature of fluorescent lamps, all shipments of Fluorescents are subject to an additional charge for over pack. We reserve the right to specify the carrier on all fluorescent shipments. Call Bulbworks for the additional freight charges which will be based on the size of your order.

*Price shown here is the Bulbworks Suggested Retail Price. The Equipment Manufacturer's List Price for the same item may vary widely from what is indicated here.

For Equipment To Lamp Information, Please Consult The Bulbworks, Inc. Cross Reference Guide

Many Other Lamp Types Are Available! Please Call For Any Lamp Not Listed!

Bulbworks, Inc. CALL: 800-334-2852 (973-584-7171) • FAX: 800-335-2858 (973-584-0300)

E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com

**Replacement Lamps For
Dental Equipment**

Pictures Are For Reference Only And Are Not To Scale



**Figure Numbers Correspond To The Bulbworks Cross Reference Guide
And Price Sheet For Dental Equipment July, 2016
No Part Of This Catalog/Price Sheet May Be Reproduced Without Permission
Bulbworks, Inc. CALL: 800-334-2852 (973-584-7171) • FAX: 800-335-2858 (973-584-0300)
E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com**



**Figure Numbers Correspond To The Bulbworks Cross Reference Guide
And Price Sheet For Dental Equipment July, 2016**

No Part Of This Catalog/Price Sheet May Be Reproduced Without Permission

Bulbworks, Inc. CALL: 800-334-2852 (973-584-7171) • FAX: 800-335-2858 (973-584-0300)

E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com

Equipment To Lamp Cross Reference Guide

July, 2016

	Manufacturer	Model	Mfr's Part No.	Description	Bulbworks Part No.▲	Net Price	Fig.
A	AcuCam-New Im.	AcuCam		12V,100W	BW.EFP	\$23.15	1
	AcuCam-New Im.	AcuCam 250		120V,250W	BW.ETJ	24.97	1
	AcuCam-New Im.	Concept IV	602457	12V,75W	BW.64450AX	10.50	15
	AcuCam-New Im.	Polo		12V,100W	BW.EVA	14.15	12
	ADC	Quala		14V,35W	BW.13165-IMP	20.75	3
	A-Dec	Excellence (Older Style)	041-117-00	3.5V,2.16W	BW.041-117-00	58.20	39
				3.5V,.75A	BW.041-317-00	53.25	39
	A-Dec	"Newer Style" Handpiece (Can Use Standard Bulb or LED)	041-317-00	LED Version – Rated 3.5V,.75A	BW.041-31700LED Can Be Used To Replace Standard 041-317-00 Bulb	66.95	48
	A-Dec	Performer	041-513-00	12V,55W	BW.JA12V-55WD	7.95	28
	A-Dec	Performer	041-501-00	9W Compact Fluor.	BW.CF9DS.E.841	17.95	60
	A-Dec–All Except Performer	6300 (Older Model); Cascade	041-179-01	24V,150W	BW.FCS	7.69	14
	A-Dec	6300 (Newer Model); 571	041-709-00	17V,95W	BW.041-709-00IM	19.75	16
	A-Dec	RA24 Roto-Quick Coupler	0.03523010		BW.0.03523010	86.75	43
	ADT	Ultra Cam III		12V,50W	BW.64602	26.45	12
	Air Techniques	Vista Cam	77121	15V,150W	BW.EFR	23.15	1
	Air Techniques	Day Loader 105A		120V,7W RED	BW.7C7-RED	5.26	25
	Algotec	Walklite		8V,50W	BW.EFM	21.85	2
	Almore	OPII	3442	12V,75W	BW.JCRM12V-75W	30.15	4
	Ampco	Duo Light		12V,50W	BW.JC12V-50WG1	13.05	12
	Astron Dental	Curing Oven XL3		12V,20W	BW.FTD	9.18	5
B	Belmont	Belmont Lamp	800-01006	12V,55W	BW.JA12V-55WD	7.95	28
	Belmont	Clesta, 501, 502, 503	501-704	24V,60W	BW.JA24V60WDBT	29.50	29
	Bisco		V10105A	12V,75W	BW.JCRM12V-75W	30.15	4
	Bosworth	Curing Oven	LCO	120V,150W	BW.JCR120V150WB	31.75	1
	Bosworth	Superlight	92102	15V,150W	BW.EFR	23.15	1
Boyd	Exam Light; OpLite S1500		12V,20W	BW.ESX***	10.25	1	
C	Castle	Aseptec 2000	516983	24V,75W	BW.JCD24V-75W	32.25	19
	Castle	Aseptec 2000	516767	24V,60W	BW.JCD24V-60WBA	31.89	19
	Castle	Older Operatory		120V,300W	BW.300SB.IF	30.72	26
	Castle	Panovision		120V,125W	BW.125T10.P	32.25	20
	Castle	Rolux II	12-1909-20	30V,50W	BW.ENZ	28.65	1
	Castle	2410C	49455	120V,150W	BW.ESP	23.40	19
	Castle	Also See MDT					
	L.D.Caulk	Max Light	644725	14V,35W	BW.13165-IMP	20.75	3
	L.D.Caulk	Prisma Lite, PR1	644101	24V,200W	BW.EJL	21.85	1
	L.D.Caulk	Prismetics Lite; EU 001	644358	12V,100W	BW.JCRM12V-100W	35.60	4
	L.D.Caulk	Prolite	644815	14V,35W	BW.13165-IMP	20.75	3
	L.D.Caulk	Spectrum; 200R	644751	10V,52W	BW.3367-5267**	52.11	3
	L.D.Caulk	See Dentsply For Additional Models					
	Chayes/Virginia	Chayes Light		12V,55W	BW.JA12V-55WD	7.95	28
	Clive Craig	Desk Top Viewer		14W Fluorescent	BW.F14T12.D	15.35	63
	Clive Craig	Wall Viewer		14W Fluorescent	BW.F14T8.D	10.30	62
	Coe Laboratories	COE-Lite	4070110	15V,150W	BW.EFR	23.15	1
	Coltene/Whaledent	Coltolux	C7906	14V,35W	BW.13165-IMP	20.75	3
	Coltene/Whaledent	Coltolux II; III; IV; 75	C7925	12V,75W	BW.JCRM12V-75W	30.15	4
	Coltene/Whaledent	DI-500 Oven	C7292	12V,50W	BW.JC12V-50WG1	13.05	12
Cuda Products	M-150; MT-150	SYC0029	21V,150W	BW.EJA	18.15	1	
Cuda Products	M2-150; M2-150T; M3-150; M3-150T; M2-150-300	SY0023	21V,150W	BW.DNF	28.65	6	
Cygnus Imaging	Cygnascope		120V,250W	BW.ETJ	24.97	1	
Cygnus Imaging	Stylus 240; 2000		12V,75W	BW.EFN	23.15	1	
D	Dabi Atlante Amer.	Focus L		12V,55W	BW.JA12V-55WD	7.95	28
	Dabi Atlante Amer.	MSL 350 Handpiece	42001-0212	Generic - 1 Pk	BW.42001021GEN1	35.85	42
	Dabi Atlante Amer.	MSL 350 Handpiece	42001-0212	Generic - 3 Pk	BW.42001021GEN3	96.75	42
	Danserau Health			12V,55W	BW.JA12V-55WD	7.95	28

▲ Multiple Brands Other Than Brand Shown May Be Available – See Price List
 **Same As Caulk/Dentsply #644751 - also fits the 3M XL1000 operated at different voltage.
 *** Some Model Boyd Exam Lights Require 13.8V Which Is No Longer Available.

Call Bulbworks • 800-334-2852 (973-584-7171) For Any Manufacturer Or Model Not Listed!
E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com

Equipment To Lamp Cross Reference Guide

July, 2016

	Manufacturer	Model	Mfr's Part No.	Description	Bulbworks Part No.▲	Net Price	Fig.	
D	Darby	Superdent 2001;2500; 8001	950-9023	12V,75W	BW.JCRM12V-75W	\$30.15	4	
	Darby	Superdent 5000; 8000	950-9021	120V,150W	BW.ESD	31.80	1	
	Darby	Superdent 7500	950-9024	21V,150W	BW.EJV	23.65	1	
	DCI	Top Mount Light	1248	24V, 100W	BW.64638-IMP	16.52	14	
	Degussa	Degulux	2352-3301	12V,75W	BW.JCRM12V-75W	30.15	4	
	Demco Grinder	E96/B1	1073	12.8V,1.80A	BW.1073	1.62	56	
	Demetron	Optilux 35; 50; 150	20919	14V,35W	BW.13165-IMP	20.75	3	
	Demetron	Optilux 101	20437	12V,75W	BW.JCRM12V-75W	30.15	4	
	Demetron	Optilux 180; 360	3367-5267	10V,52W**	BW.3367-5267**	52.11	3	
	Demetron	Optilux 400	20437	12V,75W	BW.JCRM12V-75W	30.15	4	
	Demetron	400; 401; 403 High Output	20140	12V,100W	BW.JCRM12V-100W	35.60	4	
	Demetron	Optilux 380; 400; 401; 402; 403; 500; 501	21237	12V,80W	BW.21237*	36.80	4	
	Demetron	VCL 100; 200; 300	20437	12V,75W	BW.JCRM12V-75W	30.15	4	
	Den-Mat	Ceri Light	3961	12V,100W	BW.JCRM12V-100W	35.60	4	
	Den-Mat	Marathon; Visar II	3931	21V,150W	BW.EJA	18.15	1	
	Den-Mat	Marathon II	3904	21V,150W	BW.EJV	23.65	1	
	Dental & Medical Diagnosis (DMD)	Telecam		120V,150W	BW.ESD	31.80	1	
	Dentalez	AXCS Light		24V,150W	BW.FCS	7.69	14	
	Dentalez	Colorite		12V,50W	BW.JCRM12V-50W	33.19	3	
	Dentalez	Distinction		24V,150W	BW.FDS.DZE	17.05	17	
	Dentalez	201 Desktop Viewer		14W Fluorescent	BW.F14T8.D	10.30	62	
	Dentalez	202 Wall Viewer		14W Fluorescent	BW.F14T8.D	10.30	62	
	Dental Resources	Pro Cure	P600-100	120V,150W	BW.JCR120V150WB	31.75	1	
	Dental Vision Direct	Ultra Cam		12V,50W	BW.JC12V-50W	11.34	12	
	Dentamerica	Celux Operatory		12V,55W	BW.JA12V-55WD	7.95	28	
	Dentamerica	Litex 660; 680	668	12V,75W	BW.JCRM12V-75W	30.15	4	
	Dentamerica	690		14V,35W	BW.13165-IMP	20.75	3	
	Dentech		01-11-132	24V,150W	BW.FDS.DZE	17.05	17	
	Dentech		02-50-010	24V,150W	BW.FCS	7.69	14	
	Dentech	Ultima II Fiberoptics	01-11-046	12V,100W	BW.EVA	14.15	12	
	Dentech	Headlight		4V,3.4W T3 Min.	BW.HPR53	5.78	55	
	Dentech	Wall Communicator		14V,.08A	BW.7382	3.07	37	
	Dentrix	Intraoral Camera		12V,75W	BW.EFN	23.15	1	
	Dentsply	Perspective Dental Imaging; Retrospective Treatment		12V,75W	BW.EFN	23.15	1	
	Dentsply	QHL-75	644657	12V,75W	BW.JCRM12V-75W	30.15	4	
	Dentsply	Spectrum; 200R	3367-5267	9.8V,49.9W**	BW.3367-5267**	52.11	3	
	Dentsply	Super Spectrum	13865	12V,75W	BW.13865	42.55	4	
	Dentsply	TCUIII		21V,150W	BW.EJV	23.65	1	
	Dentsply	Triad 1000; 2000	70353	120V,250W	BW.ETJ	24.97	1	
	Dentsply	Triad II; D41	95775	41V,200W	BW.95775	40.85	1	
	Dentsply	See L.D.Caulk Under The Letters "C" & "L" For Additional Models						
	Design For Vision	Durolite 3000 Headlight		120V,250W	BW.BHB	41.95	6	
	Design For Vision	Quadrallite		21V,150W	BW.ELZ	32.65	6	
	DeTrey	Prismetics		12V,75W	BW.JCRM12V-75W	30.15	4	
	Dexta	Starlite, Star II		12V,20W	BW.ESX	10.25	1	
Dolan Jenner	145		120V,150W	BW.DNE	33.12	6		
Dolan Jenner	Fiberlite 170; 176		21V,150W	BW.EJV	23.65	1		
E	EFOS	EFOS-Lite	3021	12V,100W	BW.JCRM12V-100W	35.60	4	
	EFOS	EFOS-35	3026	14V,35W	BW.13165-IMP	20.75	3	
	Elipar	See Premier Dental						
	Eilman	Conversion		120V,250W	BW.EHT	20.15	33	
	Envision	Polaris Intraoral		100V,150W	BW.ESY	28.47	13	
	ESPE	See Premier Dental & 3M						
Eurodent	Eurodent Light		24V,150W	BW.FCS	7.69	14		

▲ Multiple Brands Other Than Brand Shown May Be Available – See Price List

* Some Demetron Models Were Originally equipped With 12V, 75W Lamps. Replaced With 12V, 80W Lamps.

** Same As Caulk/Dentsply #644751 - also fits the 3M XL1000 operated at different voltage.

Call Bulbworks • 800-334-2852 (973-584-7171) For Any Manufacturer Or Model Not Listed!
E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com

Equipment To Lamp Cross Reference Guide

July, 2016

	Manufacturer	Model	Mfr's Part No.	Description	Bulbworks Part No.▲	Net Price	Fig.
F	Faro	Koolbeam		24V,150W	BW.FDS.DZE	\$17.05	17
	Faro	Original Sunlight 70		24V,150W	BW.FCS	7.69	14
	Faro	P287 Polymerizing		14V,35W	BW.13165-IMP	20.75	3
	Faro	S-90 Sunlight; S2000		24V,150W	BW.FDS.DZE	17.05	17
	Faro	Upgraded Sunlight 70		24V,150W	BW.FDS.DZE	17.05	17
	Fiad	The Light		12V,55W	BW.JA12V-55WD	7.95	28
	First Medica	Hilux 100		10V,52W	BW.3367-5267	52.11	3
	First Medica	Hilux 200;250;300;350	950-302	12V,75W	BW.JCRM12V-75W	30.15	4
	First Medica	Hilux 601		12V,80W	BW.21237	36.80	4
	Forest	9000		24V,150W	BW.FDS.DZE	17.05	17
	Forest	9020	0014-081	12V55W	BW.JA12V-55WD	7.95	28
	Forest	3110 ; 9030	0014-187	12V,50W	BW.64440	8.30	12
Fuji Optical	Denta Cam 2; 5; 11; 50; 100		15V,150W	BW.EFR	23.15	1	
Fuji Optical	Vision Plus		100V,150W	BW.ESY	28.47	13	
G	Gendex	AcuCam Concept IV	602457	12V,75W	BW.64450AX	10.50	15
	Gendex	Pan X-Ray	Alignment Light	6V,0.67A	BW.HPR40	2.75	55
	Getz	Focus Activator	4501	12V,100W	BW.JCRM12V-100W	35.60	4
	Goodlite	Halogen Headlight	HG880-1000	12V,14W	BW.789	18.35	34
	Goodlite	700 Headlight	IN447-700	6.8V,1.91A	BW.88	3.10	57
H	Healthco	Aristocrat 1 (Gun)	HCOB400491	12V,75W	BW.JCRM12V-75W	30.15	4
	Healthco	Aristocrat II (Gun)	HCOB400351	14V,35W	BW.13165-IMP	20.75	3
	Healthco	Aristocrat Light		12V,55W	BW.JA12V-55WD	7.95	28
	Healthco	Celebrity		12V,75W	BW.JCRM12V-75W	30.15	4
	Healthco	DE 100; DE 200 Viewer		14W Fluorescent	BW.F14T8.D	10.30	62
	Healthco	Visible Curing Wand		15V,150W	BW.EFR	23.15	2
	Healthco	VLC Illuminator Wand	HCOB400481	120V,250W	BW.ETJ	24.97	1
	Healthco	VLC High Intensity	HCOB400711	15V,150W	BW.EFR	23.15	1
	Henry Schein	Crusader, Starbrite	100-58-35	24V,150W	BW.FDS.DZE	17.05	17
	Henry Schein	VCL-300		12V,100W	BW.JCRM12V-100W	35.60	4
	Henry Schein	See Schein For Other Models					
	I	Hirox	KH2700R		60V,21W MR11	BW.M21E001	372.00
IDT		Oral Vision-Newer		82V,85W	BW.ESJ	36.52	1
IDT		Oral Vision-Older		120V,250W	BW.ETJ	24.97	1
Integra Medical		Viper Cam		12V,100W	BW.EFP-OSR	25.45	2
Interlight		LSM Lite Specialty Metalworks		24V,150W	BW.FCS	7.69	14
J	Ivoclar	See Vivadent					
	Jeneric Pentron	Avante GT 332* Newer	M44BA	12V,75W	BW.JCRM12V-75W	30.15	4
	Jeneric Pentron	Some Models		120V,150W	BW.JCR120V150WB	31.75	1
	J. Morita	J. Morita Light		24V,150W	BW.FCS	7.69	14
	J. Morita	Jetlite 4000 Plus	24-970302	12V,75W	BW.JCRM12V-75W	30.15	4
	J. Morita	Jetlite 4000; 69363	24-980302	10V,52W	BW.3367-5267	52.11	3
K	Johnson & Johnson	Fotofil		12V,75W	BW.EFN	23.15	1
	Kavo	Handpiece Lux Control Module	553-3881 Kavo Brand	3.3V,4W 2 Piece Pack	BW.553-3881 Kavo Brand	166.75	42
	Kavo	Handpiece Models 196; 200; 700; 701; 702TLC	1002-2928 Kavo Brand	3.3V,4W 2 Piece Pack	BW.553-3881 Kavo Brand	166.75	42
	Kavo	Handpiece Lux Control Module; Handpiece Models 196; 200; 700; 701; 702TLC (Can Use Standard Bulb or LED)	1002-2928 Kavo Brand	3.3V,4W	BW.5533881-GEN1 Generic Brand - Single	38.60	42
				3.3V,4W	BW.5533881-GEN3 Generic Brand-3 Pack	106.50	42
				LED Version - Rated 3.3V,4W	BW.553-3881-LED Can Be Used To Replace Standard 553-3881 Bulb	66.95	47
	Kavo	Optica	025-5001	12V,100W	BW.EFP	23.15	1
	Kavo	Polar II; Vicon DLS; Vicon LS; Vicon LS II (Older)	555-5001	120V,250W	BW.ETJ	24.97	1
	Kavo	Continued On The Next Page					

▲ Multiple Brands Other Than Brand Shown May Be Available – See Price List

*Older Models Used JCR12V-75W MR8 Available from the Equipment Manufacturer Only Part #M44AA (Round Pins).

Call Bulbworks • 800-334-2852 (973-584-7171) For Any Manufacturer Or Model Not Listed!

E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com

Equipment To Lamp Cross Reference Guide

July, 2016

	Manufacturer	Model	Mfr's Part No.	Description	Bulbworks Part No.▲	Net Price	Fig.
K	Kavo	Continued From The Previous Page					
	Kavo	Vicon	555-3081	6V,20W	BW.M34	\$14.90	11
	Kavo	Vicon LS	905-8500	12V,100W	BW.EVA	14.15	12
	Kavo	Vicon LS II (Newer)	905-8500	12V,100W	BW.EVA	14.15	12
	Kerr Manufacturing	Command Unit	15512	21V,150W	BW.EJM	26.65	1
	Kinetic Instruments	Sunlite II	1450	21V,150W	BW.EJV	23.65	1
	Kinetic Instruments	Sunspot	1511	12V,100W	BW.JCRM12V-100W	35.60	4
	Kinetic Instruments	Kinetic Handpiece (Can Use Standard Bulb or LED)	1222	LED Version - Rated 3.5V	BW.1222-LED Can Be Used To Replace Standard 1222 Bulb	68.25	44
	Knight Manufacturing		148656	24V,150W	BW.FDS.DZE	17.05	17
	Knight Manufacturing	Model L		24V,150W	BW.FCS	7.69	14
	Kulzer	Translux CL	2244	12V,100W	BW.JCRM12V-100W	35.60	4
	Kulzer	Translux EC; EC-S	2045	15V,150W	BW.EFR	23.15	1
	Kulzer	Translux II	2045	15V,150W	BW.EFR	23.15	1
L	L.D. Caulk	Max Light; 100	644725	14V,35W	BW.13165-IMP	20.75	3
	L.D. Caulk	Spectrum; 200R	644751	10V,52W**	BW.3367-5267**	52.11	3
	L.D. Caulk	Prisma Lite; PR1	644101	24V,200W	BW.EJL	21.85	1
	L.D. Caulk	Prismetics Lite	644358	12V,100W	BW.JCRM12V-100W	35.60	4
	L.D. Caulk	Prolite	644815	14V,35W	BW.13165-IMP	20.75	3
	L.D. Caulk	See Dentsply For Additional Models					
	Lares	Apollo	10276	12V,75W	BW.JCRM12V-75W	30.15	4
	Lares	Quartz Fire SX	10037	120V,250W	BW.ETJ	24.97	1
	Lasac Tech	LTI-5000		12V,75W	BW.JCRM12V-75W	30.15	4
	Lester A Dine			120V,150W	BW.JCD120V150WB	26.42	12
	Lester A Dine	Oral Dine		21V,150W	BW.EJA	18.15	1
	Light Tech	Headlight		6V,10W	BW.787	14.02	34
	Litema	Arcus II		10V,52W	BW.3367-5267	52.11	3
	Litema	Pluraflex HL150		15V,150W	BW.EFR	23.15	1
	LSM Lite Specialty Metalworks	Interlight		24V,150W	BW.FCS	7.69	14
	Luxo	Viewer		15W Fluorescent	BW.F15T8.D	10.60	62
	M	Marus	Patient Light		24V,150W	BW.FDS.DZE	17.05
Megadental		Megalux		12V,75W	BW.JCRM12V-75W	30.15	4
Mega Physik		Cromalux 100		12V,75W	BW.JCRM12V-75W	30.15	4
Mega Physik		Optima EC 250	EHJ	24V,250W	BW.EHJ	12.75	14
MDT		Castle Aseptec 2000	516983	24V,75W*	BW.JCD24V-75W*	32.25	19
MDT		Excalibur		12V,55W	BW.JA12V-55WD	7.95	28
MDT		Rolux I	12-2024-10	24V,200W	BW.EJL	21.85	1
MDT		Rolux II	12-1909-20	30V,50W	BW.ENZ	28.65	1
Micro Motors		LHS 710 Fiberoptic	LHS711	12V,100W	BW.EVA	14.15	12
Midmark		See Knight Manufacturing and Ritter					
Midmark		153585	002-1154-00	17V,95W	BW.041-709-00IM	11.90	16
Midmark		Castle Aseptec 2000	516983	24V,60W	BW.JCD24V-60WBA	31.89	19
Midwest		Oral Illuminator II; III	490709	21V,150W	BW.ELZ	32.65	6
Midwest		Power Optic Coupler Stylus Single Handpiece (Can Use Standard Bulb or LED)	790257	1 Piece Pack	BW.790257GEN1	63.95	36
Midwest				LED Version – Rated 3.5V	BW.790257-LED Can Be Used To Replace Standard 790257 Bulb	66.95	46
Midwest		Power Optic Handpiece	650197	3 Pack-Midwest	BW.650197	Call	35
Midwest		Power Optic In Hose Swivel	650064	2 Pack-Generic	BW.650064-GEN2	61.40	38
Midwest		Tradition I; Insight I & II	490874	120V,250W	BW.ETJ	24.97	1
Miltex		Denlite Chair Site Mirror	DP1020	1 Piece Pack	BW.DP1020	48.95***	36
Miltex	Denlite Portable Mirror	DP1520	1 Piece Pack	BW.DP1520	36.45***	36	
MK Dent	QC5016 Handpiece Coupler	BU7012	Generic 1 Pack	BW.BU7012-GEN1	38.60	42	
MK Dent	QC5016 Handpiece Coupler	BU7012	Generic 3 Pack	BW.BU7012-GEN3	106.50	42	
J.Morita	J.Morita Light		24V,150W	BW.FCS	7.69	14	
Motloid	Compolite 660		12V,75W	BW.JCRM12V-75W	30.15	4	
N	New Image	AcuCam		120V,250W	BW.ETJ	24.97	1
	New Image	Dental View		100V,150W	BW.ESY	28.47	13
	New Image	Intra View		12V,100W	BW.EFP-OSR	25.45	2
	New Image	Multi Cam	600188	12V,100W	BW.EFP-OSR	25.45	2

▲ Multiple Brands Other Than Brand Shown May Be Available – See Price List

*Original Lamp May Have Been A 24V, 60W (MDT#516767) – See BW.JCD24V-60WBA

**Also Fits The 3M XL1000 Operated At Different Voltage.

***Price Is For Stock On Hand Only

Call Bulbworks • 800-334-2852 (973-584-7171) For Any Manufacturer Or Model Not Listed!

E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com

Equipment To Lamp Cross Reference Guide

July, 2016

	Manufacturer	Model	Mfr's Part No.	Description	Bulbworks Part No.▲	Net Price	Fig.
N	NSK	High Speed Handpiece; Low Speed Handpiece; NL-10 Quick Connect Swivel Coupler	Y900529 Generic Brand	For NSK Handpiece	BW.Y900529-GEN1 Generic Brand - 1 Pk	Call	36
	NSK	High Speed Handpiece; Low Speed Handpiece; NL-10 Quick Connect Swivel Coupler	Y900529 Generic Brand	For NSK Handpiece	BW.Y900529-GEN3 Generic Brand - 3 Pk	Call	36
	NSK	Novlite		12V,75W	BW.64613-IMP	\$24.85	4
	NSK	Optica Fiber Optic	U025-410	12V,100W	BW.EFP-OSR	25.45	2
	NSK	Opticure		12V,100W	BW.EFP-OSR	25.45	2
O	Oral Vision	Oral Vision		120V,250W	BW.ETJ	24.97	1
	Orascoptic Research	New Unit	Door on Side	21V,150W	BW.EJV	23.65	1
	Orascoptic Research	Older Unit	Door on Front	21V,150W	BW.EJA	18.15	1
	O'Ryan Industries	Omega 110; 510		120V,150W	BW.DNE	33.12	6
	O'Ryan Industries	Model 750	750/BULB	12V,75W	BW.JCRM12V-75W	30.15	4
P	Panasonic	KS102		12V,100W	BW.EFP	23.15	1
	Panasonic	Monitor		21V,150W	BW.EJV	23.65	1
	Patterson Dental	Dental Cam		14.5V,90W	BW.EPX	31.65	1
	Patterson Dental	Solarc, Reveal DLS400	AL1824, 09500	24W Solarc	BW.09500-USH	409.50	8
	Patterson Dental	TCL390		14V,35W	BW.13165-IMP	20.75	3
	Patterson Dental	TCL490		10V,52W	BW.3367-5267	52.11	3
	Patterson Dental	TCL490 Plus		12V,80W	BW.21237	36.80	4
	Patterson Dental	VCL300		12V,75W	BW.JCRM12V-75W	30.15	4
	Pelton & Crane	Executive Cabinet		8W Fluorescent	BW.F8T5.D	8.60	61
	Pelton & Crane	Executive Cabinet	8044	13W Fluorescent	BW.F13T5.CW	10.30	61
	Pelton & Crane	Executive Viewer		14W Fluorescent	BW.F14T8.D	10.30	62
	Pelton & Crane	Indicator Panel Light		24V,50MA Slide Lamp	BW.552450	4.95	51
	Pelton & Crane	LF-1; LF-2		25V,150W	BW.Q150T4.CL25V	19.63	30
	Pelton & Crane	LF1; LF2 Panel Lights		28V,0.04A	BW.387	2.03	53
	Pelton & Crane	Spirit		14W Fluorescent	BW.F14T12.CW	14.98	63
	Pelton & Crane	Spirit Fiber Optic	3304909	3 Piece Pack	BW.650197	Call	35
	Penn-Med	Original Wand	301204	120V,150W	BW.ESD	31.80	1
	Penn-Med	200L Wand	301203	21V,150W	BW.EJV	23.65	1
	Penn-Med	300L; 310 (Gun)	301315	12V,75W	BW.JCRM12V-75W	30.15	4
	Phasealloy	Witelite Gun	3050	12V,75W	BW.JCRM12V-75W	30.15	4
	Phasealloy	Witelite Wand	3051	120V,150W	BW.ESD	31.80	1
	Philips	X-Ray		24V,50MA Slide Lamp	BW.552450	4.95	51
	Plan Meca	Delight		24V,150W	BW.FCS	7.69	14
	Porter	Daystar S-90B		24V,150W	BW.FDS.DZE	17.05	17
	Porter	Portable Communications System; Sensinol		24V,20MA	BW.462420	3.50	50
	Premier Dental*	Anesthetic Cartridge Warmer	3S6-120V	120V,3W	BW.3S6-120V	2.70	27
	Premier Dental*	Elipar I; Visio Beta	77160	13V,75W	BW.70201102004	177.50	2
	Premier Dental*	Elipar II; Elipar Mini Spot High Light		12V,75W	BW.64613-IMP	24.85	4
	Pro Den Systems	Executor	20-055	12V,100W	BW.JCRM12V-100W	35.60	4
	Pro Den Systems	Polylite 500; 1000; 2000	20-016	120V,150W	BW.ESD	31.80	1
	Pro Den Systems	Pro-Lite DLS	72-190	12V,50W	BW.ENL	19.37	1
	Pro Den Systems	ProScope 1000	42-9035	14V,35W	BW.13165-IMP	20.75	3
Pro Den Systems	ProScope 2000;4000		12V,75W	BW.JCRM12V-75W	30.15	4	
Pro Den Tec	Oral Video Scope		120V,250W	BW.ETJ	24.97	1	
Pro Den Tec	Prism		12V,100W	BW.EFP-OSR	25.45	2	
Pro Den Tec	Pro Cam C1000		30V,50W	BW.ENZ	28.65	1	
Pro Vision Imaging	DIS; Pro Vision		120V,250W	BW.ETJ	24.97	1	
Q	Quality Aspirators	Fiber Optic Daylight		21V,150W	BW.EJA	18.15	1

▲ Multiple Brands Other Than Brand Shown May Be Available – See Price List

Call Bulbworks • 800-334-2852 (973-584-7171) For Any Manufacturer Or Model Not Listed!
E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com

Equipment To Lamp Cross Reference Guide

July, 2016

	Manufacturer	Model	Mfr's Part No.	Description	Bulbworks Part No.▲	Net Price	Fig.
R	Rembrandt	Virtuoso Phase II		24V,250W	BW.ELC	\$18.23	2
	Rinn Duplicator	Slimline Duplicator	003203	8W 12" Black Light Blue	BW.F8T5.BLB	18.25	64
	Rinn Duplicator	Viewlight Safelight		120V,7W RED	BW.7C7-RED	5.26	25
	Rinn Duplicator	Viewlight Safelight		120V,7W ORANGE	BW.7C7-ORANGE	5.26	25
	Rinn	Universal Viewboxes 67-0400;67-0440;67-1000;67-1020;67-6000;67-6020;67-8001		8W Fluorescent	BW.F8T5.D	8.60	61
	Ritter	355 Procedure Light		24V,100W	BW.64638-IMP	16.52	14
	Ritter	Model H		120V,125W	BW.125T10.P	32.25	20
	Ritter	Midmark 151		120V,25W	BW.25A-120V	3.27	21
	Ritter	Midmark 354	002-060-0500	21V,150W	BW.EKE	16.26	1
	Ritter	Starlite "J"		28V,75W	BW.Q75CL.DC-28V	20.10	19
	Ritter	Super Starlite "K"		24V,150W	BW.FDS.DZE	17.05	17
	Robin	UL1000		24V,150W	BW.FDS.DZE	17.05	17
	Rolux	MDT	12-2024-10	24V,200W	BW.EJL	21.85	1
	S	Sable	Access L & Rotomax L Handpiece	Sable 1	1 Piece Pack	BW.SABLE-1	Call
Sable		Access L & Rotomax L Handpiece	Sable 5	5 Piece Pack	BW.SABLE-5	Call	36
Satalec		Servolite		8V,50W	BW.EFM	21.85	2
Schein		Crusader; Star Brite		24V,150W	BW.FDS.DZE	17.05	17
Schein		Economy; VCL200	101-0045	12V,75W	BW.JCRM12V-75W	30.15	4
Schein		Maxima 250	101-6062	14V,35W	BW.13165-IMP	20.75	3
Schein		VLC 300	100-5657	12V,100W	BW.JCRM12V-100W	35.60	4
Schultz Dental		Elcelux		12V,100W	BW.JCRM12V-100W	35.60	4
Sci Can		Dual Lux II Curing		120V,250W	BW.ETJ	24.97	1
Sci Can		Dual Lux II; III Illum.		21V,150W	BW.EKE	16.26	1
Sci Can		Dual Lux III Curing		21V,150W	BW.EJA	18.15	1
SDS-Summit Dental		Overhead Light		24V,70W	BW.JA24V-70WD	14.90	28
Seltzer		Potato Chair	150W TLC	95V,150W	BW.MHR150N	315.00	10
Siemens		A1/E1 Sirona		24V,150W	BW.FCS	7.69	14
Siemens		Sirona Operatory		12V,50W	BW.64440	8.30	12
Siemens		X-Ray		12V,20W	BW.64115	15.10	55
Simplified Systems		Communicator		28V,.04A	BW.85	2.95	40
Solarc		Patterson Dental	AL1824, WA09500, 09500	24W Solarc	BW.09500-USH	409.50	8
Solid State Systems		Initiator	6220002-290	15V,150W	BW.EFR	23.15	1
Sony		20/20 Carry Cam		21V,150W	BW.EJA	18.15	1
S.S. White		Profile VCS	17244	12V,100W	BW.JCRM12V-100W	35.60	4
Spectraview		Intra & Extraoral Camera		12V,75W	BW.EFN	23.15	1
Springhealth Penn-Med		Original Wand		120V,150W	BW.ESD	31.80	1
Springhealth Penn-Med		Powerlite 75; 100; 310 Gun	BU75	12V,75W	BW.JCRM12V-75W	30.15	4
Springhealth Penn-Med		200L Wand		21V,150W	BW.EJV	23.65	1
Star Dental		Handpiece for Quick Connect Swivel	260986	4.5-5.5V 2 Piece Pack	BW.60986	90.65	32
Star Dental		Starlite	61351	21V,150W	BW.DNF	28.65	6
Star X-Ray	DE 94; 98; 100; 120; 150; 160; 200; SW100 Viewbox		14W Fluorescent	BW.F14T8.D	10.30	62	
Star X-Ray	DE 99		15W Fluorescent	BW.F15T8.D	10.60	62	
Star X-Ray	Endodontic Viewer		8W Fluorescent	BW.F8T5.D	8.60	61	
Star X-Ray	Mini Viewer		6W Fluorescent	BW.F6T5.D	8.25	61	
Star X-Ray	Starcure 350 EH		120V,150W	BW.JCR120V150WB	31.75	1	
Star X-Ray	TR 130 Cephalometric		8W Fluorescent	BW.F8T5.D	8.60	61	
Star X-Ray	612 Duplicator		120V,7W Orange	BW.7C7-ORANGE	5.26	25	
Steris	Examiner 10		120V,150W	BW.ESP	23.40	19	
Summit Dental Sys.	Overhead Light		24V,70W	BW.JA24V-70WD	14.90	28	

▲ Multiple Brands Other Than Brand Shown May Be Available – See Price List
* Also See 3M

Call Bulbworks • 800-334-2852 (973-584-7171) For Any Manufacturer Or Model Not Listed!
E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com

Equipment To Lamp Cross Reference Guide

July, 2016

	Manufacturer	Model	Mfr's Part No.	Description	Bulbworks Part No.▲	Net Price	Fig.
T	Teledyne Getz	Focus Activator	4501	12V,100W	BW.JCRM12V-100W	\$35.60	4
	Theta Communications	Five-Lite – Square Lenses 1972 and before		6.3V,,15A	BW.40	1.62	-
	Theta Communications	Five-Lite – Square Lenses 1973 and after		6.3V,,15A	BW.47	1.62	54
	Theta Communications	Five-Lite – Square Lenses 1997 and before		6.3V,,20A	BW.381	3.40	53
	Theta Communications	Five-Lite – Square Lenses 1998 and after		6.3V,,15A	BW.47	1.62	54
	Theta Communications	Ten-Lite – Square Lenses 1997 and before		28V,,06A	BW.376	2.82	53
	Theta Communications	Ten-Lite – Square Lenses 1998 and after		28V,,04A	BW.1819	3.65	54
	Theta Communications	Ten-Lite – Round Lenses		28V,,04A	BW.1819	3.65	54
	Theta Communications	Status System – Round Lenses of Units with 600 Part Numbers		6.3V,,15A	BW.47	1.62	54
	Theta Communications	Status System – Round Lenses of Units with 700 Part Numbers		28V,,04A	BW.1819	3.65	54
	3M	Visilux I; II	78-8045-4571-9	12V,75W	BW.JCRM12V-75W	30.15	4
	3M	Visilux 5510		12V,50W	BW.JC12V-50W	11.34	12
	3M	XL1500	3367-5267	10V,52W	BW.3367-5267*	52.11	3
	3M	XL2500; 3000	78-8045-4571-9	12V,75W	BW.JCRM12V-75W	30.15	4
	3M-ESPE	Elipar I; Elipar Visio Beta	70-2011-0200-4	13V,75W	BW.70201102004	177.50	2
3M-ESPE	Elipar II; Elipar Mini Spot High Light	70-2011-0199-8	12V,75W	BW.64613-IMP	24.85	4	
3M-ESPE	2500	78-8131-1662-7	12V,75W	BW.64613-IMP	24.85	4	
Trojan Research	Ultraeye		21V,150W	BW.EKE	16.26	1	
Trophy	Intra Oral		12V,75W	BW.EFN	23.15	1	
U	Unitek/3M	Ortholux; XT	78-8045-4571-9	12V,75W	BW.JCRM12V-75W	30.15	4
V	Vacudent			24V,150W	BW.FDS.DZE	17.05	17
	Vector	VKS6 Handpiece	VRB Bulb	Generic 1 Pack	BW.VRB-GEN1	Call	41
	Vector	VKS6 Handpiece	VRB Bulb	Generic 3 Pack	BW.VRB-GEN3	Call	41
	Vidicam			120V,150W	BW.JCR120V150WB	31.75	1
	Vivadent	Astralis V, VII		12V,75W	BW.14552	49.55	3
	Vivadent	Astralis 10	561398	13V,100W	BW.561398	77.95	3
	Vivadent	Heliolux; Heliolux II	7412, 133762	14V,35W	BW.13165-IMP	20.75	3
	Vivadent	Heliolux DLX; GT	7512	12V,75W	BW.JCRM12V-75W	30.15	4
	Vivadent	Heliomat B Single	7210	21V,150W	BW.EJA	18.15	1
	Vivadent	Heliomat H10 Single; Heliomat Multi	7110	15V,150W	BW.EFR	23.15	1
W	W & H	W & H Handpiece (Can Use Standard Bulb or LED)	22076	LED Version – Rated as 3.5V	BW.22076-LED Can Be Used To Replace Standard 22076 Bulb	66.95	45
	Weber	Panavision		120V,125W	BW.125T10.P	32.25	20
	Weber	SK81G		12V,100W	BW.JA12V-100WD	11.40	28
	Weber	SK100 Skylight		12V,100W	BW.FDT	14.73	18
	Weber	SK800; SK900G Skylight	SK580	14.5V,90W	BW.EPV	24.85	1
	Weil Dental	Elcelux		14V,35W	BW.13165-IMP	20.75	3
	Welch Allyn	DA-L-400	WA09800	21W	BW.WA09800	Call	8
	Welch Allyn	Denlite Chair Site Mirror	DP1020	1 Piece Pack	BW.DP1020	48.95**	36
	Welch Allyn	Denlite Portable Mirror	DP1520	1 Piece Pack	BW.DP1520	36.45**	36
	Welch Allyn	Exam Light II	WA04200	12V,35W	BW.EPN	38.15	1
	Welch Allyn	Solarc, Reveal DLS400	AL1824, WA09500, 09500	24W Solarc	BW.09500-USH	409.50	8
	S.S.White	Profile VCS	17244	12V,100W	BW.JCRM12V-100W	35.60	4
	Wykle Research	Witelite 1500	3052	12V,75W	BW.JCRM12V-75W	30.15	4

▲ Multiple Brands Other Than Brand Shown May Be Available – See Price List

* Same as Caulk/Dentsply #644751 operated at different voltage.

** Price Is For Stock On Hand Only

Call Bulbworks • 800-334-2852 (973-584-7171) For Any Manufacturer Or Model Not Listed!
E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com

Replacement Lamps For Dental Equipment Distributor/Industry Number Cross Reference Guide

July, 2016

Distrib. No. (Efos)	Distrib. No. (Efos & Caulk)	Distrib. No. (Chapman-Huffman)	Distrib. No. (Quala)	Philips No.	Osram No.	GE, Industry Std., Or ANSI No.	Bulbworks No.
3008	640 30080	16-432-00	-	-	-	EZK	BW.JCR120V150WB
3009	644751	16-433-00	QUA-230	13298	-	-	BW.3367-5267
3010	640 30100	-	-	-	64613	64613	BW.64613
3011	640 30110	16-360-00	QUA-160	5995	-	EJM	BW.EJM
3014	640 30140	16-330-00	QUA-165	6423	64634	EFR	BW.EFR
3015	640 30150	16-410-00	QUA-170	-	-	ETJ	BW.ETJ
3016	640 30160	16-420-00	QUA-175	13865	64617	JCRM12V-75W	BW.JCRM12V-75W
3017	640 30170	16-350-00	QUA-180	13164	64644	EJL	BW.EJL
3018	640 30180	16-340-00	QUA-185	-	-	EJA	BW.EJA
3019	640 30190	16-300-00	-	-	-	DNE	BW.DNE
3020	640 30200	16-320-00	QUA-190	6834	64627	EFP	BW.EFP
3021	640 30210	16-430-00	QUA-195	-	64624	JCRM12V-100W	BW.JCRM12V-100W
3022	640 30220	16-370-00	QUA-200	-	-	EJV	BW.EJV
3023	640 30230	-	-	-	-	DNF	BW.DNF
3024	640 30240	16-390-00	-	-	-	ELZ	BW.ELZ
3025	640 30250	-	-	-	-	-	BW.650197
3026	644725	16-431-00	QUA-205	13165	-	JCRM14V-35W	BW.13165
3027	640 30270	16-380-00	QUA-210	-	-	EKE	BW.EKE
3028	640 30280	16-400-00	QUA-215	-	-	ESD	BW.ESD
3031	640 30310	16-054-00	-	-	-	041-317-00	BW.041-317-00
3032	640 30320	-	-	-	64440	JC12V-50W	BW.JC12V-50W
3033	640 30330	-	-	-	-	-	***
3034	640 30340	16-445-00	-	-	-	ESY	BW.ESY
3035	640 30350	16-436-00	QUA-220	6853	64615	EFN	BW.EFN
3036	-	-	-	-	-	Q50MR16/FL/1 - ENL	BW.ENL
3038	640 30380	16-437-00	-	7748	64655	EHJ	BW.EHJ
3040	640 30400	16-435-00	-	6847	64607	EFM	BW.EFM
3047	640 30470	-	-	-	-	FC8T9/CW	BW.FC8T9.CW
3048	640 30480	-	-	-	-	F6T5/D	BW.F6T5.D
3049	640 30490	-	QUA-140	-	-	F8T5/D	BW.F8T5.D
3049B	640 3049B	-	-	-	-	F8T5.BLB	BW.F8T5.BLB
3050	640 30500	-	-	-	-	F6T5/CW	BW.F6T5.CW
3051	640 30510	-	-	-	-	F8T5/CW	BW.F8T5.CW
3052	640 30520	-	-	-	-	F13T5/CW	BW.F13T5.CW
3053	640 30530	-	QUA-155	-	-	F14T12/CW	BW.F14T12.CW

*** Not Available Through Bulbworks As Of April, 2004

Continued On The Next Page...

For Additional Cross Reference Information Contact
Bulbworks, Inc. • P.O. Box 586, Succasunna, NJ 07876
CALL: 800-334-2852 (973-584-7171) • FAX: 800-335-2858 (973-584-0300)
E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com

Replacement Lamps For Dental Equipment Distributor/Industry Number Cross Reference Guide

July, 2016

Continued From The Previous Page

Distrib. No. (Efos)	Distrib. No. (Efos & Caulk)	Distrib. No. (Chapman-Huffman)	Distrib. No. (Quala)	Philips No.	Osram No.	GE, Industry Std., Or ANSI No.	Bulbworks No.
3054	640 30540	-	QUA-150	-	-	F14T12/D	BW.F14T12.D
3055	640 30550	-	-	-	-	F14T8/CW	BW.F14T8.CW
3056	640 30560	-	QUA-145	-	-	F14T8/D	BW.F14T8.D
3057	640 30570	-	-	-	-	F15T8/CW	BW.F15T8.CW
3058	640 30580	-	-	-	-	F15T8/D	BW.F15T8.D
3059	640 30590	16-411-00	-	7724	64623	EVA	BW.EVA
4070	640 40700	16-020-00	QUA-115	7158	64640	JC24V-150W	BW.FCS
4071	640 40710	16-040-00	QUA-100	12336	64151	JA12V-55WD (H3)	BW.JA12V-55WD
4072	640 40720	16-035-00	QUA-105	125T10/P	-	125T10/P	BW.125T10.P
4073	640 40730	16-015-00	-	-	-	Q250/MC/CL-120V EHT	BW.EHT
4074	640 40740	16-025-00	QUA-110	5974	64643	FDS/DZE	BW.FDS.DZE
4075	640 40700	16-020-00	QUA-115	7158	64640	JC24V-150W	BW.FCS
4076	640 40760	16-050-00	QUA-120	-	-	Q150T4/CL-25V	BW.Q150T4.CL25V
4077	640 40770	16-050-00	QUA-120	-	-	Q150T4/CL-25V - GE Brand	BW.Q150T4CL25GE
4078	640 40780	-	-	-	-	-	BW.PC TOOL
4079	640 40790	16-045-00	QUA-125	-	-	Q75CL/DC-28V	BW.Q75CL.DC-28V
4080	640 40800	16-010-00	-	-	-	EPV	BW.EPV
4081	640 40810	16-030-00	-	5973	64628	FDT	BW.FDT
4082	-	-	-	-	-	100A21/SP	BW.100A21SP
4084	640 40840	16-351-00	-	-	-	ENZ	BW.ENZ
4085	640 40850	16-011-00	-	6642	-	ESX	BW.ESX
5616	-	16-320-00	QUA-190	6834	64627	EFP	BW.EFP
8023	-	-	-	-	-	FC8T9/CW	BW.FC8T9.CW
96000	640 96000	-	-	-	CF9DS/BLUE	CF9DS/BLUE	BW.CF9DS.BLUE
-	-	16-012-00	-	-	-	JCD24V-75W	BW.JCD24V-75W
-	-	16-058-00	-	-	-	650064	BW.650064-GEN2
-	-	16-061-00	QUA-135	-	-	553-3881	BW.553-3881
-	-	16-062-00	-	-	-	Y900529	BW.Y900529-GEN1
-	-	16-063-00	-	-	-	42001021	BW.42001021GEN1
-	-	16-440-00	QUA-225	-	-	21237	BW.21237
-	-	-	QUA-130	-	64451	-	BW.64451
-	-	-	QUA-131	14623	-	-	BW.041-709-00

For Additional Cross Reference Information, Contact:

Bulbworks, Inc. • P.O. Box 586, Succasunna, NJ 07876
CALL: 800-334-2852 (973-584-7171) • FAX: 800-335-2858 (973-584-0300)
E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com



Dental Waste Recycling Products



BW.SUPPLY-062 - 2 GALLON DENTAL WASTE RECYCLING KIT - \$103.95

Kit includes: 2 pails with locking lids - 1 pail (½ gallon) for up to 2 lbs. of wet dental amalgam waste and 1 pail (2 gallons) for up to 24 lbs. of other dental waste including separator filters, traps and lead foil. Also includes instructions, terms and conditions, generator certification form, return shipping box, prepaid return shipping label, proof of purchase with serial number, and online certificates of recycling.

Note: The ½ gallon pail must be placed inside the 2 gallon pail and then placed in the shipping box prior to return shipping. Bilingual – English and Spanish Packaging.

DOT Limits the amount of Elemental Mercury in the container to 1 lb. Not available in AK, HI, PR, VI, AR, RI. Not available for export to Canada.

For Additional Information Go To: <https://www.bulbworks.com/light-bulbs/SUPPLY-062>



BW.SUPPLY-063 - 2 GALLON DENTAL FIXER RECYCLING KIT - \$103.95

Kit includes: 1 pail with locking lid, 6 mil poly liner with liner tie for up to 2 gallons of dental fixer solution. Also includes instructions, terms and conditions, return shipping box, prepaid return shipping label, proof of purchase with serial number, and online certificates of recycling.

Note: The 2 gallon pail must be placed in the shipping box prior to return shipping. Bilingual – English and Spanish Packaging.

DOT Limits the amount of Elemental Mercury in the container to 1 lb. Not available in AK, HI, PR, VI, AR, RI. Not available for export to Canada.

For Additional Information Go To: <https://www.bulbworks.com/light-bulbs/SUPPLY-063>



BW.SUPPLY-125 - 5 GALLON DENTAL WASTE RECYCLING KIT - \$143.50

Kit includes: 2 pails with locking lids - 1 pail (½ gallon) for up to 2 lbs. of wet dental amalgam waste and 1 pail (5 gallons) for up to 64 lbs. of other dental waste including separator filters, traps and lead foil. Also includes instructions, terms and conditions, generator certification form, prepaid return shipping label, proof of purchase with serial number, and online certificates of recycling.

Note: The ½ gallon pail must be placed inside the 5 gallon pail prior to return shipping. Bilingual – English and Spanish Packaging.

DOT Limits the amount of Elemental Mercury in the container to 1 lb. Not available in AK, HI, PR, VI, AR, RI. Not available for export to Canada.

For Additional Information Go To: <https://www.bulbworks.com/light-bulbs/SUPPLY-125>



BW.SUPPLY-234 – 11.3 GALLON SEPARATOR FILTER PAIL - \$187.50

Includes: 1 pail with locking lid, 6 mil poly liner with liner tie for up to 70 lbs. dental filters and chair trap waste. Also includes instructions, terms and conditions, prepaid return shipping label, proof of purchase with serial number, and online certificates of recycling.

Note: Bilingual – English and Spanish Packaging.

DOT Limits the amount of Elemental Mercury in the container to 1 lb. Not available in AK, HI, PR, VI, AR, RI. Not available for export to Canada.

For Additional Information Go To: <https://www.bulbworks.com/light-bulbs/SUPPLY-234>



Bulbworks, Inc. - P.O. Box 586 - Succasunna, NJ 07876

Phone: 800-334-2852

Fax: 800-335-2858

Email: bulbwork@bulbworks.com

www.bulbworks.com

Getting To Know Us... (Terms & Conditions)

ORDERING IS EASY...

Call Us TOLL FREE • 800-334-BULB(2852) • Our business hours are 8:30 AM to 5:00 PM Eastern Time Monday through Thursday and 8:30 AM to 4:00 PM on Friday. Outside of business hours? Use our after hours message machine and your call will be returned during the next business day or you may...

Fax Us ON OUR 24 HOUR TOLL FREE FAX LINE • 800-335-2858 • All orders received by fax are acknowledged promptly either by phone or return fax.

Order Online • www.bulbworks.com • All orders placed online will be acknowledged promptly by return e-mail.

No Minimum Order • We appreciate your business no matter how small or large your order may be. Orders with a merchandise value of under \$10.00 will be subject to a Small Order Handling Charge of \$3.00 (In addition to Merchandise and Shipping).

Fast Shipment • Orders received by 3:00 PM Eastern Time will be **shipped the same day** that we receive them. Orders are normally shipped by UPS but all express services are available. Bulbworks carefully inspects all lamps prior to shipment and we use the best packaging methods available.

Terms • We accept MasterCard, Visa, Discover, & American Express. Our terms for *approved* open accounts are 2% 10 Days, Net 30, freight prepaid and added to the invoice. We would be glad to open an account for your convenience. Just ask!



Order With Confidence • All items carry a full warranty. Bulbworks sells only top quality brand name products - the same brands used by the equipment manufacturers. Should product failure, breakage, or defect occur, contact us for immediate replacement!

Original Manufacturer's Packaging • Each lamp is shipped in its original packaging from the lamp manufacturer – Eiko, Philips, GE, Ushio, Osram - Sylvania, Welch Allyn, and others. Part numbers are clearly labeled and are easy to read to eliminate confusion at reorder time. Helpful equipment reference information appears on every carton.

All Prices Are Subject To Change Without Notice. Some Prices May Be For Stock On Hand Only.
No Part Of This Catalog May Be Reproduced Without Permission.



P.O. Box 586 Succasunna, NJ 07876
CALL: 800-334-BULB (2852) • (973-584-7171)
FAX: 800-335-2858 • (973-584-0300)

E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com



**Need Something Besides Dental Bulbs? Stop Looking & Call Us!
We Have What You Need! We Stock All Types Of
Bulbs, Sockets, and Batteries Including:**

- **MICROSCOPE & SCIENTIFIC BULBS** Replacement Lamps For American Optical, Bausch & Lomb, Jena, Leitz, Nikon, Olympus, Swift, Unitron, Zeiss, And Others. Lamps For Biological, Biomedical, Educational, Forensic, Industrial, Metallurgical, Research, And Other Applications. Fluorescent And LED Microscope Illuminators.
- **MEDICAL BULBS** Replacement Lamps For Operatory, Surgical, Clinical, Fiber Optic, Laser, Analytical, Veterinary, X-Ray Viewing, And Other Applications. Lamps For All Types Of Medical Equipment.
- **INSPECTION** Replacement Lamps For Optical Comparators Including Baty Gagemaster, J & L, Micro Vu, Mitutoyo, Nikon, OGP, Scherr Tumico, Vermont Precision, And Others. Replacement Lamps for Microscopes And Other Inspection Devices.
- **GRAPHIC ARTS BULBS** Replacement Lamps For AB Dick, Agfa, Amergraph, Brown, Burgess, Douthitt, Itek, Kenro, NuArc, Olec, Sandmar, Screen Theimer, 3M, And More.
- **OPTICAL AND OPHTHALMIC BULBS** Replacement Lamps For American Optical, Bausch & Lomb, Copeland, Haag Streit, Marco, Mentor, Nikon, Reichert, Topcon, Welch Allyn, And Others. Lamps For Optical And Ophthalmic Equipment Including Lensometers, Ophthalmoscopes, Projectors, Retinoscopes, Slit Lamps, Tonometers, And Others.
- **GENERAL AND SPECIALTY LIGHTING** Replacement Lamps And LEDs For Color Enhancement And Correction, Commercial, Decorative, Energy Saving, General Service, Long Life, Office, Problem Solving, Rough Service, And Other Applications.
- **GERMICIDAL, BLACK LIGHT, AND OTHER UV BULBS** Replacement Lamps For Air And Water Purification. For Use In Hospitals, Laboratories, Manufacturing, Food Handling, And Many Other Applications.
- **PROJECTOR AND AUDIO VISUAL BULBS** Replacement Lamps For Overhead Projectors, Slide And Movie Projectors, Microfilm and Microfiche Readers, Micrographic And Other Applications.
- **SOCKETS AND LAMP HOLDERS** Replacement Sockets For All Lamp Types Including Halogen, Incandescent, Vacuum, & Others. Sockets From Bender & Wirth, Buhl, Gilway, Christiana Industries, Cole Hersee, Gilway, Sylvania, Ushio, And Others.
- **BATTERIES** Batteries For All Types Of Devices. Alkaline, Lithium, Nickel Metal Hydride, Heavy Duty, Zinc Air, Rechargeable, Industrial, Hearing Aid, Photo, Electronic, Medical, Keyless Entry, Watch and Others. All Sizes Available.

We Supply All Specialty Lamp Types Including: Fluorescent, Incandescent, Halogen, Mercury Vapor, Metal Halide, Miniature, Quartz, Xenon, And Other Specialty Bulbs.

BULBWORKS Supplies Only Top Quality Brand Names - The Same Brands Used By The Equipment Manufacturers. Lamps From Eiko, GE, Hosobuchi, Narva, Nikon, Opti Quip, Osram Sylvania, Philips, Ushio, Welch Allyn, ORC, Voltarc, Zeiss, And Many Other Specialty Lamp Manufacturers

Call BULBWORKS for Your Next Replacement Lamp Requirement!

Bulbworks, Inc. • P.O. Box 586 • Succasunna, NJ 07876

CALL: 800-334-BULB(2852) • (973-584-7171)

FAX: 800-335-2858 • (973-584-0300)

E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com