



## LED And Fluorescent Microscope Illuminators



**BW.HX-1**  
(Fluorescent)

110V, 8W, Daylight, 6400 Degrees K, Ring Shape, Fluorescent Ring Light Assembly For Microscopes, Self Contained Unit Includes Bulb, Power Supply And Connecting Cable, Maximum Connecting Caliber Is 60mm, Includes 1 Additional Bulb  
Uses Replacement Bulb BW.92ADL

**Net Price: \$69.50**



**BW.HX-2**  
(Fluorescent)

The Fluorescent Version Of This Item BW.HX-2 Has Been Discontinued. The Nearest Replacement Is The LED Version BW.HXD-7. Replacement Fluorescent Bulbs For HX-2 Remain Available As Our Part Number BW.FCB92D.

**Discontinued**  
See BW.HXD-7



**BW.HX-4**  
(Fluorescent)

120V, 8W, Daylight, 6400 Degrees K, Ring Shape, Fluorescent Ring Light For Microscopes, Self Contained Unit Includes Bulb, Power Supply, Switch And Connecting Cable, Maximum Connecting Caliber Is 55mm  
Uses Replacement Bulb BW.FCB80A

**Net Price: \$69.50**



**BW.HXD-1**  
(LED)

120V, 6W, Daylight, 6000 Degrees K, Ring Shape, LED Ring Light For Microscopes, 60 LEDs, 96-265 VAC, Dimmable, Includes Plastic Power Supply With Controller, Maximum Connecting Caliber Is 60mm, Lamp Head Is ABS Plastic And Aluminum

**Net Price: \$89.00**



**BW.HXD-2**  
(LED)

120V, 12W, 3200 Degrees K, Gooseneck Microscope Illuminator With Self Contained Power Supply, 2 Tubes, LED, 6W Per Tube, 300mm (11.8") Gooseneck Length, 92-220 VAC

**Net Price: \$109.50**



**BW.HXD-2B**  
(LED)

120V, 16W, Daylight, 6000 Degrees K, Gooseneck Microscope Illuminator With Self Contained Power Supply And Heavy Duty Metal Base, 500mm (19.7") Gooseneck Length, Dual Head, LED, Each Head 8 Watts, Dimmable, 50-60 HZ, 85-265 VAC

**Net Price: \$159.50**



**BW.HXD-3**  
(LED)

120V, 6W, Daylight, 6000 Degrees K, Ring Shape, Ring Light For Microscope, Includes Ring Light, Heavy Duty Metal Power Supply With Controller And Connector, 60 LEDs (6 Watt), Dimmable, 96-220 VAC, Maximum Connecting Caliber Is 60mm, Lamp Head Is Sturdy Plastic And Metal

**Net Price: \$116.95**



**BW.HXD-4**  
(LED)

120V, 3W, Daylight, 6500-7000 Degrees K, Ring Light For Microscope, Self Contained Unit includes Power Supply with Line Cord, Switch and Ring Light, Dimmable, Maximum Connecting Caliber is 60mm, Durable Plastic Housing

**Net Price: \$36.95**



**BW.HXD-5**  
(LED)

120V, 6W, Daylight, 6000 Degrees K, Ring Shape, Ring Light For Microscope, 60 LEDs (6 Watt), 96-220 VAC, Dimmable, Maximum Connecting Caliber Is 60mm, Plastic Power Supply Case, Lamp Head Is Durable Plastic And Aluminum

**Net Price: \$96.50**



**BW.HXD-6**  
(LED)

12W, Daylight, 6000 Degrees K, Ring Shape, Heavy Duty Ring Light For Microscope, Multi Level Lighting With Up To 12 Watts (80 LEDs), 5 Levels Of Lighting Intensity - Full or Additional Sub Area Lighting 1/2, 1/4, 1/8, Max. Connecting Caliber Is 60mm, Heavy Duty Power Supply & Lamp Head

**Net Price: \$179.50**



**BW.HXD-7**  
(LED)

120V, 6W, Daylight, 6000 Degrees K, Ring Shape, Ring Light For Microscope, Self Contained Unit Includes Power Supply, Switch and Ring Light, Dimmable, Maximum Connecting Caliber Is 60mm, Durable Plastic Housing

**Net Price: \$99.50**



**BW.HXD-9**  
(LED)

120V, 16W, Gooseneck Microscope Illuminator For Mounting Directly On Scope, Including Mounting Brackets, Dual LED Light Heads Each 3-8 Watts, Dimmable, 590mm (23.2") Gooseneck Length, 50-60 HZ, 85-265 VAC

**Net Price: Call**



**BW.HXP-295**  
(LED)

120V, 4W, Circular Shape, LED Stage Illuminator For Microscope, 95mm Diameter x 4.88mm Thick, 4 Watt Assembly Includes Power Supply And Connector, Light Is For Bottom Lighting, Provides Transmitted Light, Dimmable, 90-265 VAC

**Net Price: \$69.00**



**BW.92ADL**  
(Fluorescent)

8 Watt Replacement Lamp for HX-1 Microscope Illuminator

**Net Price: \$59.50**



**BW.FCB92D**  
(Fluorescent)

10 Watt Replacement Lamp for HX-2 Microscope Illuminator

**Net Price: \$115.40**



**BW.FCB80A**  
(Fluorescent)

8 Watt Replacement Lamp for HX-4 Microscope Illuminator

**Net Price: \$44.75**