

BULBWORKS

SPECIALTY LIGHTBULBS AND SOCKETS

Cross Reference Guide For Replacement Lamps For DENTAL EQUIPMENT Consumer



Replacement Lamps For Operatory, Curing, Fiberoptic,
Handpiece, Examination, Video, X-ray Viewing & More

Bulbworks, Inc.
P.O. Box 586, Succasunna, NJ 07876
Call: **800-334-2852** · (973-584-7171)
Fax: **800-335-2858** · (973-584-0300)
E-mail: bulbwork@bulbworks.com
Visit Us Online: www.bulbworks.com



Cross Reference Guide Replacement Lamps For Dental Equipment

Replacement Lamps For Operatory, Curing, Fiberoptic,
Handpiece, Examination, Video, X-Ray Viewing & More

Welcome...

| | |
|--|---------|
| Pictures (Figure Table)..... | Page 1 |
| Equipment To Lamp Cross Reference Guide..... | Page 3 |
| Distributor/Industry Number Cross Reference Guide..... | Page 10 |
| Terms & Conditions..... | Page 12 |

Pay Less, Expect More...

Just because our lamps cost less doesn't mean you'll get less! Learn why Bulbworks has earned the business of Dental Dealers and Equipment Manufacturers all over the country. We will quickly become a valuable asset to *your* business.

Bulbworks offers only the highest quality specialty lamps for use in Dental Equipment. We feature lamps from Ushio, Philips, Osram/Sylvania, GE, Eiko, and many other top manufacturers - the same brands used by the equipment manufacturers, but you'll save big!

You can buy with confidence - All lamps carry a full warranty!
We will immediately replace any defective lamp!

We Go Out Of Our Way For You...

At Bulbworks, customer service is not a department - *it's our philosophy!* Expect knowledgeable, courteous, and memorable service. For reliability, we can't be beat! With Toll Free Phone and Fax Numbers, we're convenient, economical, and *we ship fast!*

Bulbworks, Inc. • P.O. Box 586 • Succasunna, NJ 07876
CALL: 800-334-2852 • (973-584-7171) • TOLL-FREE 24-HOUR FAX: 800-335-2858 (973-584-0300)
E-mail: bulbwork@bulbworks.com • Visit Us Online: www.bulbworks.com



Replacement Lamps For Dental Equipment

Pictures Are For Reference Only And Are Not To Scale



Figure 1



Figure 2



Figure 3



Figure 4



Figure 5




Figure 6



Figure 7



Figure 8



Figure 9



Figure 10



Figure 11



Figure 12



Figure 13



Figure 14



Figure 15



Figure 16



Figure 17



Figure 18



Figure 19



Figure 20



Figure 21



Figure 22



Figure 23



Figure 24



Figure 25



Figure 26



Figure 27



Figure 28



Figure 29



Figure 30



Figure 31



Figure 32



Figure 33

Figure Numbers Correspond To The Bulbworks Cross Reference Guide
No Part Of This Catalog May Be Reproduced Without Permission
Bulbworks, Inc. CALL: 800-334-2852 (973-584-7171) • FAX: 800-335-2858 (973-584-0300)
E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com



Figure 34



Figure 35



Figure 36



Figure 37



Figure 38



Figure 39



Figure 40



Figure 41



Figure 42



Figure 43



Figure 44



Figure 45



Figure 46



Figure 47



Figure 48



Figure 49



Figure 50



Figure 51



Figure 52



Figure 53



Figure 54



Figure 55



Figure 56



Figure 57



Figure 58



Figure 59



Figure 60



Figure 61

Figure 62

Figure 63



Figure 64



Figure 65

Figure Numbers Correspond To The Bulbworks Cross Reference Guide
No Part Of This Catalog/Price Sheet May Be Reproduced Without Permission
Bulbworks, Inc. CALL: 800-334-2852 (973-584-7171) • FAX: 800-335-2858 (973-584-0300)
E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com

Equipment To Lamp Cross Reference Guide
For Pricing See Separate Bulbworks Price Sheet

| | Manufacturer | Model | Mfr's Part No. | Description | Bulbworks Part No.▲ | Fig. | |
|----------------|----------------------------|--|-----------------------|-------------------------------|--|----------------|----|
| A | AcuCam-New Im. | AcuCam | | 12V,100W | BW.EFP | 1 | |
| | AcuCam-New Im. | AcuCam 250 | | 120V,250W | BW.ETJ | 1 | |
| | AcuCam-New Im. | Concept IV | 602457 | 12V,75W | BW.64450AX | 15 | |
| | AcuCam-New Im. | Polo | | 12V,100W | BW.EVA | 12 | |
| | | Quala | | 14V,35W | BW.13165-IMP | 3 | |
| | A-Dec | Excellence (Older Style) | 041-117-00 | 3.5V,2.16W | BW.041-117-00 | 39 | |
| | | | | 3.5V,.75A | BW.041-317-00 | 39 | |
| | A-Dec | "Newer Style" Handpiece (Can Use Standard Bulb or LED) | 041-317-00 | LED Version – Rated 3.5V,.75A | BW.041-31700LED Can Be Used To Replace Standard 041-317-00 Bulb | 48 | |
| | A-Dec | Performer | 041-513-00 | 12V,55W | BW.JA12V-55WD | 28 | |
| | A-Dec | Performer | 041-501-00 | 9W Compact Fluor. | BW.CF9DS.E.841 | 60 | |
| | A-Dec–All Except Performer | 6300 (Older Model); Cascade | 041-179-01 | 24V,150W | BW.FCS | 14 | |
| | A-Dec | 6300 (Newer Model); 571 | 041-709-00 | 17V,95W | BW.041-709-00IM | 16 | |
| | A-Dec | RA24 Roto-Quick Coupler | 0.03523010 | | BW.0.03523010 | 43 | |
| | ADT | Ultra Cam III | | 12V,50W | BW.64602 | 12 | |
| | Air Techniques | Vista Cam | 77121 | 15V,150W | BW.EFR | 1 | |
| | Air Techniques | Day Loader 105A | | 120V,7W RED | BW.7C7-RED | 25 | |
| | Algotec | Walklite | | 8V,50W | BW.EFM | 2 | |
| | Almore | OPII | 3442 | 12V,75W | BW.JCRM12V-75W | 4 | |
| | Ampco | Duo Light | | 12V,50W | BW.JC12V-50WG1 | 12 | |
| | Astron Dental | XL2, XL3, XLGT | | 9W Blue | BW.CF9DS.BLUE | 59 | |
| | Astron Dental | Curing Oven XL3 | | 12V,20W | BW.FTD | 5 | |
| | B | Belmont | Belmont Lamp | 800-01006 | 12V,55W | BW.JA12V-55WD | 28 |
| | | Belmont | Clesta, 501, 502, 503 | 501-704 | 24V,60W | BW.JA24V60WDBT | 29 |
| Bisco | | | V10105A | 12V,75W | BW.JCRM12V-75W | 4 | |
| Bosworth | | Curing Oven | LCO | 120V,150W | BW.JCR120V150WB | 1 | |
| Bosworth | | Superlight | 92102 | 15V,150W | BW.EFR | 1 | |
| Boyd | | Exam Light; OpLite S1500 | | 12V,20W | BW.ESX*** | 1 | |
| C | Castle | Aseptic 2000 | 516983 | 24V,75W | BW.JCD24V-75W | 19 | |
| | Castle | Aseptic 2000 | 516767 | 24V,60W | BW.JCD24V-60WBA | 19 | |
| | Castle | Older Operator | | 120V,300W | BW.300SB.IF | 26 | |
| | Castle | Panovision | | 120V,125W | BW.125T10.P | 20 | |
| | Castle | Rolux II | 12-1909-20 | 30V,50W | BW.ENZ | 1 | |
| | Castle | 2410C | 49455 | 120V,150W | BW.ESP | 19 | |
| | Castle | Also See MDT | | | | | |
| | L.D.Caulk | Max Light | 644725 | 14V,35W | BW.13165-IMP | 3 | |
| | L.D.Caulk | Prisma Lite, PR1 | 644101 | 24V,200W | BW.EJL | 1 | |
| | L.D.Caulk | Prismetics Lite; EU 001 | 644358 | 12V,100W | BW.JCRM12V-100W | 4 | |
| | L.D.Caulk | Prolite | 644815 | 14V,35W | BW.13165-IMP | 3 | |
| | L.D.Caulk | Spectrum; 200R | 644751 | 10V,52W | BW.3367-5267** | 3 | |
| | L.D.Caulk | See Dentsply For Additional Models | | | | | |
| | Chayes/Virginia | Chayes Light | | 12V,55W | BW.JA12V-55WD | 28 | |
| | Clive Craig | Desk Top Viewer | | 14W Fluorescent | BW.F14T12.D | 63 | |
| | Clive Craig | Wall Viewer | | 14W Fluorescent | BW.F14T8.D | 62 | |
| | Coe Laboratories | COE-Lite | 4070110 | 15V,150W | BW.EFR | 1 | |
| | Coltene/Whaledent | Coltolux | C7906 | 14V,35W | BW.13165-IMP | 3 | |
| | Coltene/Whaledent | Coltolux II; III; IV; 75 | C7925 | 12V,75W | BW.JCRM12V-75W | 4 | |
| | Coltene/Whaledent | DI-500 Oven | C7292 | 12V,50W | BW.JC12V-50WG1 | 12 | |
| | Cuda Products | M-150; MT-150 | SYC0029 | 21V,150W | BW.EJA | 1 | |
| | Cuda Products | M2-150; M2-150T; M3-150; M3-150T; M2-150-300 | SY0023 | 21V,150W | BW.DNF | 6 | |
| | Cygnus Imaging | Cygnascope | | 120V,250W | BW.ETJ | 1 | |
| Cygnus Imaging | Stylus 240; 2000 | | 12V,75W | BW.EFN | 1 | | |
| D | Dabi Atlante Amer. | Focus L | | 12V,55W | BW.JA12V-55WD | 28 | |
| | Dabi Atlante Amer. | MSL 350 Handpiece | 42001-0212 | Generic - 1 Pk | BW.42001021GEN1 | 42 | |
| | Dabi Atlante Amer. | MSL 350 Handpiece | 42001-0212 | Generic - 3 Pk | BW.42001021GEN3 | 42 | |
| | Danserau Health | | | 12V,55W | BW.JA12V-55WD | 28 | |

▲ Multiple Brands Other Than Brand Shown May Be Available – See Price List

**Same As Caulk/Dentsply #644751 - also fits the 3M XL1000 operated at different voltage.

*** Some Model Boyd Exam Lights Require 13.8V Which Is No Longer Available.

Call Bulbworks • 800-334-2852 (973-584-7171) For Any Manufacturer Or Model Not Listed!

E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com

EUPS2020A DENT

Equipment To Lamp Cross Reference Guide
For Pricing See Separate Bulbworks Price Sheet

| | Manufacturer | Model | Mfr's Part No. | Description | Bulbworks Part No.▲ | Fig. | |
|----------|----------------------------------|---|----------------|-----------------|---------------------|-----------------|---|
| D | Darby | Superdent 2001;2500; 8001 | 950-9023 | 12V,75W | BW.JCRM12V-75W | 4 | |
| | Darby | Superdent 5000; 8000 | 950-9021 | 120V,150W | BW.ESD | 1 | |
| | Darby | Superdent 7500 | 950-9024 | 21V,150W | BW.EJV | 1 | |
| | DCI | Top Mount Light | 1248 | 24V, 100W | BW.64638-IMP | 14 | |
| | Degussa | Degulux | 2352-3301 | 12V,75W | BW.JCRM12V-75W | 4 | |
| | Demco Grinder | E96/B1 | 1073 | 12.8V,1.80A | BW.1073 | 56 | |
| | Demetron | Optilux 35; 50; 150 | 20919 | 14V,35W | BW.13165-IMP | 3 | |
| | Demetron | Optilux 101 | 20437 | 12V,75W | BW.JCRM12V-75W | 4 | |
| | Demetron | Optilux 180; 360 | 3367-5267 | 10V,52W** | BW.3367-5267** | 3 | |
| | Demetron | Optilux 400 | 20437 | 12V,75W | BW.JCRM12V-75W | 4 | |
| | Demetron | 400; 401; 403 High Output | 20140 | 12V,100W | BW.JCRM12V-100W | 4 | |
| | Demetron | Optilux 380; 400; 401; 402; 403; 500; 501 | 21237 | 12V,80W | BW.21237* | 4 | |
| | Demetron | VCL 100; 200; 300 | 20437 | 12V,75W | BW.JCRM12V-75W | 4 | |
| | Den-Mat | Ceri Light | 3961 | 12V,100W | BW.JCRM12V-100W | 4 | |
| | Den-Mat | Marathon; Visar II | 3931 | 21V,150W | BW.EJA | 1 | |
| | Den-Mat | Marathon II | 3904 | 21V,150W | BW.EJV | 1 | |
| | Dental & Medical Diagnosis (DMD) | Telecam | | 120V,150W | BW.ESD | 1 | |
| | Dentalez | AXCS Light | | 24V,150W | BW.FCS | 14 | |
| | Dentalez | Colorite | | 12V,50W | BW.JCRM12V-50W | 3 | |
| | Dentalez | Distinction | | 24V,150W | BW.FDS.DZE | 17 | |
| | Dentalez | 201 Desktop Viewer | | 14W Fluorescent | BW.F14T8.D | 62 | |
| | Dentalez | 202 Wall Viewer | | 14W Fluorescent | BW.F14T8.D | 62 | |
| | Dental Resources | Pro Cure | P600-100 | 120V,150W | BW.JCR120V150WB | 1 | |
| | Dental Vision Direct | Ultra Cam | | 12V,50W | BW.JC12V-50W | 12 | |
| | Dentamerica | Celux Operatory | | 12V,55W | BW.JA12V-55WD | 28 | |
| | Dentamerica | Litex 660; 680 | 668 | 12V,75W | BW.JCRM12V-75W | 4 | |
| | Dentamerica | Luxus 550; 580 | | 9W Blue | BW.CF9DS.BLUE | 59 | |
| | Dentamerica | 690 | | 14V,35W | BW.13165-IMP | 3 | |
| | Dentech | | 01-11-132 | 24V,150W | BW.FDS.DZE | 17 | |
| | Dentech | | 02-50-010 | 24V,150W | BW.FCS | 14 | |
| | Dentech | Ultima II Fiberoptics | 01-11-046 | 12V,100W | BW.EVA | 12 | |
| | Dentech | Headlight | | 4V,3.4W T3 Min. | BW.HPR53 | 55 | |
| | Dentech | Wall Communicator | | 14V,.08A | BW.7382 | 37 | |
| | Dentrix | Intraoral Camera | | 12V,75W | BW.EFN | 1 | |
| | Dentsply | Perspective Dental Imaging; Retrospective Treatment | | 12V,75W | BW.EFN | 1 | |
| | Dentsply | QHL-75 | 644657 | 12V,75W | BW.JCRM12V-75W | 4 | |
| | Dentsply | Spectrum; 200R | 3367-5267 | 9.8V,49.9W** | BW.3367-5267** | 3 | |
| | Dentsply | Super Spectrum | 13865 | 12V,75W | BW.13865 | 4 | |
| | Dentsply | TCUIII | | 21V,150W | BW.EJV | 1 | |
| | Dentsply | Triad 1000; 2000 | 70353 | 120V,250W | BW.ETJ | 1 | |
| | Dentsply | Triad II; D41 | 95775 | 41V,200W | BW.95775 | 1 | |
| | Dentsply | See L.D.Caulk Under the Letters "C" & "L" For Additional Models | | | | | |
| | Design For Vision | Durolite 3000 Headlight | | 120V,250W | BW.BHB | 6 | |
| | Design For Vision | Quadralite | | 21V,150W | BW.ELZ | 6 | |
| | DeTrey | Prismetics | | 12V,75W | BW.JCRM12V-75W | 4 | |
| | Dexta | Starlite, Star II | | 12V,20W | BW.ESX | 1 | |
| | Dolan Jenner | 145 | | 120V,150W | BW.DNE | 6 | |
| | Dolan Jenner | Fiberlite 170; 176 | | 21V,150W | BW.EJV | 1 | |
| | E | EFOS | EFOS-Lite | 3021 | 12V,100W | BW.JCRM12V-100W | 4 |
| | | EFOS | EFOS-35 | 3026 | 14V,35W | BW.13165-IMP | 3 |
| Elipar | | See Premier Dental | | | | | |
| Ellman | | Conversion | | 120V,250W | BW.EHT | 33 | |
| Envision | | Polaris Intraoral | | 100V,150W | BW.ESY | 13 | |
| ESPE | | See Premier Dental & 3M | | | | | |
| Eurodent | | Eurodent Light | | 24V,150W | BW.FCS | 14 | |

▲ Multiple Brands Other Than Brand Shown May Be Available – See Price List

* Some Demetron Models Were Originally equipped With 12V, 75W Lamps. Replaced With 12V, 80W Lamps.

** Same As Caulk/Dentsply #644751 - also fits the 3M XL1000 operated at different voltage.

Call Bulbworks • 800-334-2852 (973-584-7171) For Any Manufacturer Or Model Not Listed!

E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com

Equipment To Lamp Cross Reference Guide
For Pricing See Separate Bulbworks Price Sheet

| | Manufacturer | Model | Mfr's Part No. | Description | Bulbworks Part No.▲ | Fig. |
|-----------------|-------------------|--|----------------------|-----------------------------|--|--------|
| F | Faro | Koolbeam | | 24V,150W | BW.FDS.DZE | 17 |
| | Faro | Original Sunlight 70 | | 24V,150W | BW.FCS | 14 |
| | Faro | P287 Polymerizing | | 14V,35W | BW.13165-IMP | 3 |
| | Faro | S-90 Sunlight; S2000 | | 24V,150W | BW.FDS.DZE | 17 |
| | Faro | Upgraded Sunlight 70 | | 24V,150W | BW.FDS.DZE | 17 |
| | Fiad | The Light | | 12V,55W | BW.JA12V-55WD | 28 |
| | First Medica | Hilux 100 | | 10V,52W | BW.3367-5267 | 3 |
| | First Medica | Hilux 200;250;300;350 | 950-302 | 12V,75W | BW.JCRM12V-75W | 4 |
| | First Medica | Hilux 601 | | 12V,80W | BW.21237 | 4 |
| | Forest | 9000 | | 24V,150W | BW.FDS.DZE | 17 |
| | Forest | 9020 | 0014-081 | 12V55W | BW.JA12V-55WD | 28 |
| | Forest | 3110 ; 9030 | 0014-187 | 12V,50W | BW.64440 | 12 |
| | Fuji Optical | Denta Cam 2; 5; 11; 50; 100 | | 15V,150W | BW.EFR | 1 |
| Fuji Optical | Vision Plus | | 100V,150W | BW.ESY | 13 | |
| G | Gendex | AcuCam Concept IV | 602457 | 12V,75W | BW.64450AX | 15 |
| | Gendex | Pan X-Ray | Alignment Light | 6V,0.67A | BW.HPR40 | 55 |
| | Getz | Focus Activator | 4501 | 12V,100W | BW.JCRM12V-100W | 4 |
| | Goodlite | Halogen Headlight | HG880-1000 | 12V,14W | BW.789 | 34 |
| | Goodlite | 700 Headlight | IN447-700 | 6.8V,1.91A | BW.88 | 57 |
| H | Healthco | Aristocrat 1 (Gun) | HCOB400491 | 12V,75W | BW.JCRM12V-75W | 4 |
| | Healthco | Aristocrat II (Gun) | HCOB400351 | 14V,35W | BW.13165-IMP | 3 |
| | Healthco | Aristocrat Light | | 12V,55W | BW.JA12V-55WD | 28 |
| | Healthco | Celebrity | | 12V,75W | BW.JCRM12V-75W | 4 |
| | Healthco | DE 100; DE 200 Viewer | | 14W Fluorescent | BW.F14T8.D | 62 |
| | Healthco | Visible Curing Wand | | 15V,150W | BW.EFR | 2 |
| | Healthco | VLC Illuminator Wand | HCOB400481 | 120V,250W | BW.ETJ | 1 |
| | Healthco | VLC High Intensity | HCOB400711 | 15V,150W | BW.EFR | 1 |
| | Henry Schein | Crusader, Starbrite | 100-58-35 | 24V,150W | BW.FDS.DZE | 17 |
| | Henry Schein | VCL-300 | | 12V,100W | BW.JCRM12V-100W | 4 |
| | Henry Schein | See Schein For Other Models | | | | |
| | Hirox | KH2700R | | 60V,21W MR11 | BW.M21E001 | 9 |
| | I | IDT | Oral Vision-Newer | | 82V,85W | BW.ESJ |
| IDT | | Oral Vision-Older | | 120V,250W | BW.ETJ | 1 |
| Integra Medical | | Viper Cam | | 12V,100W | BW.EFP-OSR | 2 |
| Interlight | | LSM Lite Specialty Metalworks | | 24V,150W | BW.FCS | 14 |
| Ivoclar | | See Vivadent | | | | |
| J | Jeneric Pentron | Avante GT332 Newer | M44BA | 12V,75W | BW.JCRM12V-75W | 4 |
| | Jeneric Pentron | Avante GT332 Older | M44AA | 12V,75W | BW.JCR12V75WMR8 | - |
| | Jeneric Pentron | Some Models | | 120V,150W | BW.JCR120V150WB | 1 |
| | J. Morita | J. Morita Light | | 24V,150W | BW.FCS | 14 |
| | J. Morita | Jetlite 4000 Plus | 24-970302 | 12V,75W | BW.JCRM12V-75W | 4 |
| | J. Morita | Jetlite 4000; 69363 | 24-980302 | 10V,52W | BW.3367-5267 | 3 |
| | Johnson & Johnson | Fotofil | | 12V,75W | BW.EFN | 1 |
| K | Kavo | Handpiece Lux Control Module | 553-3881 Kavo Brand | 3.3V,4W 2 Piece Pack | BW.553-3881 Kavo Brand | 42 |
| | Kavo | Handpiece Models 196; 200; 700; 701; 702TLC | 1002-2928 Kavo Brand | 3.3V,4W 2 Piece Pack | BW.553-3881 Kavo Brand | 42 |
| | | | | 3.3V,4W | BW.5533881-GEN1 Generic Brand - Single | 42 |
| | | | | 3.3V,4W | BW.5533881-GEN3 Generic Brand-3 Pack | 42 |
| | Kavo | Handpiece Lux Control Module; Handpiece Models 196; 200; 700; 701; 702TLC (Can Use Standard Bulb or LED) | 1002-2928 Kavo Brand | LED Version - Rated 3.3V,4W | BW.553-3881-LED Can Be Used To Replace Standard 553-3881 Bulb | 47 |
| | Kavo | Optica | 025-5001 | 12V,100W | BW.EFP | 1 |
| | Kavo | Polar II; Vicon DLS; Vicon LS; Vicon LS II (Older) | 555-5001 | 120V,250W | BW.ETJ | 1 |
| | Kavo | Continued On The Next Page | | | | |

▲ Multiple Brands Other Than Brand Shown May Be Available – See Price List

Call Bulbworks • 800-334-2852 (973-584-7171) For Any Manufacturer Or Model Not Listed!

E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com

Equipment To Lamp Cross Reference Guide
For Pricing See Separate Bulbworks Price Sheet

| | Manufacturer | Model | Mfr's Part No. | Description | Bulbworks Part No.▲ | Fig. |
|--------------------------------------|-----------------------------|---|-----------------|---|--|---------------------|
| K | Kavo | Continued From Previous Page | | | | |
| | Kavo | Vicon | 555-3081 | 6V,20W | BW.JC6V-20W | 11 |
| | Kavo | Vicon LS | 905-8500 | 12V,100W | BW.EVA | 12 |
| | Kavo | Vicon LS II (Newer) | 905-8500 | 12V,100W | BW.EVA | 12 |
| | Kerr Manufacturing | Command Unit | 15512 | 21V,150W | BW.EJM | 1 |
| | Keystone | Dental Curing | | 9W Black Light | BW.DS9W.78 | 59 |
| | Kinetic Instruments | Sunlite II | 1450 | 21V,150W | BW.EJV | 1 |
| | Kinetic Instruments | Sunspot | 1511 | 12V,100W | BW.JCRM12V-100W | 4 |
| | Kinetic Instruments | Kinetic Handpiece (Can Use Standard Bulb or LED) | 1222 | LED Version - Rated 3.5V | BW.1222-LED Can Be Used To Replace Standard 1222 Bulb | 44 |
| | Knight Manufacturing | | 148656 | 24V,150W | BW.FDS.DZE | 17 |
| | Knight Manufacturing | Model L | | 24V,150W | BW.FCS | 14 |
| | Kulzer | Translux CL | 2244 | 12V,100W | BW.JCRM12V-100W | 4 |
| | Kulzer | Translux EC; EC-S | 2045 | 15V,150W | BW.EFR | 1 |
| | Kulzer | Translux II | 2045 | 15V,150W | BW.EFR | 1 |
| | L | L.D. Caulk | Max Light; 100 | 644725 | 14V,35W | BW.13165-IMP |
| L.D. Caulk | | Spectrum; 200R | 644751 | 10V,52W** | BW.3367-5267** | 3 |
| L.D. Caulk | | Prisma Lite; PR1 | 644101 | 24V,200W | BW.EJL | 1 |
| L.D. Caulk | | Prismetics Lite | 644358 | 12V,100W | BW.JCRM12V-100W | 4 |
| L.D. Caulk | | Prolite | 644815 | 14V,35W | BW.13165-IMP | 3 |
| L.D. Caulk | | See Dentsply For Additional Models | | | | |
| Lares | | Apollo | 10276 | 12V,75W | BW.JCRM12V-75W | 4 |
| Lares | | Quartz Fire SX | 10037 | 120V,250W | BW.ETJ | 1 |
| Lasac Tech | | LTI-5000 | | 12V,75W | BW.JCRM12V-75W | 4 |
| Lester A Dine | | | | 120V,150W | BW.JCD120V150WB | 12 |
| Lester A Dine | | Oral Dine | | 21V,150W | BW.EJA | 1 |
| Light Tech | | Headlight | | 6V,10W | BW.787 | 34 |
| Litema | | Arcus II | | 10V,52W | BW.3367-5267 | 3 |
| Litema | | Pluraflex HL150 | | 15V,150W | BW.EFR | 1 |
| LSM Lite Specialty Metalworks | | Interlight | | 24V,150W | BW.FCS | 14 |
| Luxo | Viewer | | 15W Fluorescent | BW.F15T8.D | 62 | |
| M | Marus | Patient Light | | 24V,150W | BW.FDS.DZE | 17 |
| | Megadental | Megalux | | 12V,75W | BW.JCRM12V-75W | 4 |
| | Mega Physik | Cromalux 100 | | 12V,75W | BW.JCRM12V-75W | 4 |
| | Mega Physik | Optima EC 250 | EJH | 24V,250W | BW.EJH | 14 |
| | MDT | Castle Aseptec 2000 | 516983 | 24V,75W* | BW.JCD24V-75W* | 19 |
| | MDT | Excalibur | | 12V,55W | BW.JA12V-55WD | 28 |
| | MDT | Rolux I | 12-2024-10 | 24V,200W | BW.EJL | 1 |
| | MDT | Rolux II | 12-1909-20 | 30V,50W | BW.ENZ | 1 |
| | Micro Motors | LHS 710 Fiberoptic | LHS711 | 12V,100W | BW.EVA | 12 |
| | Midmark | 153585 | 002-1154-00 | 17V,95W | BW.041-709-00IM | 16 |
| | Midmark | Castle Aseptec 2000 | 516983 | 24V,60W | BW.JCD24V-60WBA | 19 |
| | Midmark | See Knight Manufacturing and Ritter | | | | |
| | Midwest | Oral Illuminator II; III | 490709 | 21V,150W | BW.ELZ | 6 |
| | Midwest | Power Optic Coupler Stylus Single Handpiece (Can Use Standard Bulb or LED) | 790257 | 1 Piece Pack LED Version - Rated 3.5V | BW.790257GEN1 BW.790257-LED Can Be Used To Replace Standard 790257 Bulb | 36 46 |
| | Midwest | Power Optic Handpiece | 650197 | 3 Pack-Midwest | BW.650197 | 35 |
| Midwest | Power Optic In Hose Swivel | 650064 | 2 Pack-Generic | BW.650064-GEN2 | 38 | |
| Midwest | Tradition I; Insight I & II | 490874 | 120V,250W | BW.ETJ | 1 | |
| Miltex | Denlite Chair Site Mirror | DP1020 | 1 Piece Pack | BW.DP1020 | 36 | |
| Miltex | Denlite Portable Mirror | DP1520 | 1 Piece Pack | BW.DP1520 | 36 | |
| MK Dent | QC5016 Handpiece Coupler | BU7012 | Generic 1 Pack | BW.BU7012-GEN1 | 42 | |
| MK Dent | QC5016 Handpiece Coupler | BU7012 | Generic 3 Pack | BW.BU7012-GEN3 | 42 | |
| J.Morita | J.Morita Light | | 24V,150W | BW.FCS | 14 | |
| Motloid | Compolite 660 | | 12V,75W | BW.JCRM12V-75W | 4 | |
| N | New Image | AcuCam | | 120V,250W | BW.ETJ | 1 |
| | New Image | Dental View | | 100V,150W | BW.ESY | 13 |
| | New Image | Intra View | | 12V,100W | BW.EFP-OSR | 2 |
| | New Image | Multi Cam | 600188 | 12V,100W | BW.EFP-OSR | 2 |

▲ Multiple Brands Other Than Brand Shown May Be Available – See Price List
*Original Lamp May Have Been A 24V, 60W (MDT#516767) – See BW.JCD24V-60WBA
**Also Fits The 3M XL1000 Operated At Different Voltage.

Equipment To Lamp Cross Reference Guide
For Pricing See Separate Bulbworks Price Sheet

| | Manufacturer | Model | Mfr's Part No. | Description | Bulbworks Part No.▲ | Fig. |
|---------------------------|--|---|--------------------------|------------------------|--|------|
| N | NSK | High Speed Handpiece; Low Speed Handpiece; NL-10 Quick Connect Swivel Coupler | Y900529 Generic Brand | For NSK Handpiece | BW.Y900529-GEN1 Generic Brand - 1 Pk | 36 |
| | NSK | High Speed Handpiece; Low Speed Handpiece; NL-10 Quick Connect Swivel Coupler | Y900529 Generic Brand | For NSK Handpiece | BW.Y900529-GEN3 Generic Brand - 3 Pk | 36 |
| | NSK | Novlite | | 12V,75W | BW.64613-IMP | 4 |
| | NSK | Optica Fiber Optic | U025-410 | 12V,100W | BW.EFP-OSR | 2 |
| | NSK | Opticure | | 12V,100W | BW.EFP-OSR | 2 |
| O | Oral Vision | Oral Vision | | 120V,250W | BW.ETJ | 1 |
| | Orascoptic Research | New Unit | Door on Side | 21V,150W | BW.EJV | 1 |
| | Orascoptic Research | Older Unit | Door on Front | 21V,150W | BW.EJA | 1 |
| | O'Ryan Industries | Omega 110; 510 | | 120V,150W | BW.DNE | 6 |
| | O'Ryan Industries | Model 750 | 750/BULB | 12V,75W | BW.JCRM12V-75W | 4 |
| | Panasonic | KS102 | | 12V,100W | BW.EFP | 1 |
| Panasonic | Monitor | | 21V,150W | BW.EJV | 1 | |
| Patterson Dental | Dental Cam | | 14.5V,90W | BW.EPX | 1 | |
| Patterson Dental | Solarc, Reveal DLS400 | AL1824, 09500 | 24W Solarc | BW.09500-USH | 8 | |
| Patterson Dental | TCL390 | | 14V,35W | BW.13165-IMP | 3 | |
| Patterson Dental | TCL490 | | 10V,52W | BW.3367-5267 | 3 | |
| Patterson Dental | TCL490 Plus | | 12V,80W | BW.21237 | 4 | |
| Patterson Dental | VCL300 | | 12V,75W | BW.JCRM12V-75W | 4 | |
| Pelton & Crane | Executive Cabinet | | 8W Fluorescent | BW.F8T5.D | 61 | |
| Pelton & Crane | Executive Cabinet | 8044 | 13W Fluorescent | BW.F13T5.CW | 61 | |
| Pelton & Crane | Executive Viewer | | 14W Fluorescent | BW.F14T8.D | 62 | |
| Pelton & Crane | Indicator Panel Light | | 24V,50MA Slide Lamp | BW.552450 | 51 | |
| Pelton & Crane | LF-1; LF-2 | | 25V,150W | BW.Q150T4.CL25V | 30 | |
| Pelton & Crane | LF1; LF2 Panel Lights | | 28V,0.04A | BW.387 | 53 | |
| Pelton & Crane | Spirit | | 14W Fluorescent | BW.F14T12.CW | 63 | |
| Pelton & Crane | Spirit Fiber Optic | 3304909 | 3 Piece Pack | BW.650197 | 35 | |
| Penn-Med | Original Wand | 301204 | 120V,150W | BW.ESD | 1 | |
| Penn-Med | 200L Wand | 301203 | 21V,150W | BW.EJV | 1 | |
| Penn-Med | 300L; 310 (Gun) | 301315 | 12V,75W | BW.JCRM12V-75W | 4 | |
| Phasealloy | Witelite Gun | 3050 | 12V,75W | BW.JCRM12V-75W | 4 | |
| Phasealloy | Witelite Wand | 3051 | 120V,150W | BW.ESD | 1 | |
| Philips | X-Ray | | 24V,50MA Slide Lamp | BW.552450 | 51 | |
| Plan Meca | Delight | | 24V,150W | BW.FCS | 14 | |
| Porter | Daystar S-90B | | 24V,150W | BW.FDS.DZE | 17 | |
| Porter | Portable Communications System; Sensinol | | 24V,20MA | BW.462420 | 50 | |
| Premier Dental* | Anesthetic Cartridge Warmer | 3S6-120V | 120V,3W | BW.3S6-120V | 27 | |
| Premier Dental* | Elipar I; Visio Beta | 77160 | 13V,75W | BW.70201102004 | 2 | |
| Premier Dental* | Elipar II; Elipar Mini Spot High Light | | 12V,75W | BW.64613-IMP | 4 | |
| Pro Den Systems | Executor | 20-055 | 12V,100W | BW.JCRM12V-100W | 4 | |
| Pro Den Systems | Polylite 500; 1000; 2000 | 20-016 | 120V,150W | BW.ESD | 1 | |
| Pro Den Systems | Pro-Lite DLS | 72-190 | 12V,50W | BW.ENL | 1 | |
| Pro Den Systems | ProScope 1000 | 42-9035 | 14V,35W | BW.13165-IMP | 3 | |
| Pro Den Systems | ProScope 2000;4000 | | 12V,75W | BW.JCRM12V-75W | 4 | |
| Pro Den Tec | Oral Video Scope | | 120V,250W | BW.ETJ | 1 | |
| Pro Den Tec | Prism | | 12V,100W | BW.EFP-OSR | 2 | |
| Pro Den Tec | Pro Cam C1000 | | 30V,50W | BW.ENZ | 1 | |
| Pro Vision Imaging | DIS; Pro Vision | | 120V,250W | BW.ETJ | 1 | |
| Q | Quality Aspirators | Fiber Optic Daylight | | 21V,150W | BW.EJA | 1 |

▲ Multiple Brands Other Than Brand Shown May Be Available – See Price List

*Also See 3M

Call Bulbworks • 800-334-2852 (973-584-7171) For Any Manufacturer Or Model Not Listed!

E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com

EUPS2020A DENT

Equipment To Lamp Cross Reference Guide
For Pricing See Separate Bulbworks Price Sheet

| | Manufacturer | Model | Mfr's Part No. | Description | Bulbworks Part No.▲ | Fig. |
|-----------------------|----------------------|---|--------------------------------|--------------------------|---------------------|------------|
| R | Rembrandt | Virtuoso Phase II | | 24V,250W | BW.ELC | 2 |
| | Rinn Duplicator | Slimline Duplicator | 003203 | 8W 12" Black Light Blue | BW.F8T5.BLB | 64 |
| | Rinn Duplicator | Viewlight Safelight | | 120V,7W RED | BW.7C7-RED | 25 |
| | Rinn Duplicator | Viewlight Safelight | | 120V,7W ORANGE | BW.7C7-ORANGE | 25 |
| | Rinn | Universal Viewboxes 67-0400;67-0440;67-1000;67-1020;67-6000;67-6020;67-8001 | | 8W Fluorescent | BW.F8T5.D | 61 |
| | Ritter | 355 Procedure Light | | 24V,100W | BW.64638-IMP | 14 |
| | Ritter | Model H | | 120V,125W | BW.125T10.P | 20 |
| | Ritter | Midmark 151 | | 120V,25W | BW.25A-120V | 21 |
| | Ritter | Midmark 354 | 002-060-0500 | 21V,150W | BW.EKE | 1 |
| | Ritter | Starlite "J" | | 28V,75W | BW.Q75CL.DC-28V | 19 |
| | Ritter | Super Starlite "K" | | 24V,150W | BW.FDS.DZE | 17 |
| | Robin | UL1000 | | 24V,150W | BW.FDS.DZE | 17 |
| | Rolux | MDT | 12-2024-10 | 24V,200W | BW.EJL | 1 |
| | S | Sable | Access L & Rotomax L Handpiece | Sable 1 | 1 Piece Pack | BW.SABLE-1 |
| Sable | | Access L & Rotomax L Handpiece | Sable 5 | 5 Piece Pack | BW.SABLE-5 | 36 |
| Satalec | | Servolite | | 8V,50W | BW.EFM | 2 |
| Schein | | Crusader; Star Brite | | 24V,150W | BW.FDS.DZE | 17 |
| Schein | | Economy; VCL200 | 101-0045 | 12V,75W | BW.JCRM12V-75W | 4 |
| Schein | | Maxima 250 | 101-6062 | 14V,35W | BW.13165-IMP | 3 |
| Schein | | VLC 300 | 100-5657 | 12V,100W | BW.JCRM12V-100W | 4 |
| Schultz Dental | | Elcelux | | 12V,100W | BW.JCRM12V-100W | 4 |
| Sci Can | | Dual Lux II Curing | | 120V,250W | BW.ETJ | 1 |
| Sci Can | | Dual Lux II; III Illum. | | 21V,150W | BW.EKE | 1 |
| Sci Can | | Dual Lux III Curing | | 21V,150W | BW.EJA | 1 |
| SDS-Summit Dental | | Overhead Light | | 24V,70W | BW.JA24V-70WD | 28 |
| Seltzer | | Potelt Chair | 150W TLC | 95V,150W | BW.MHR150N | 10 |
| Siemens | | A1/E1 Sirona | | 24V,150W | BW.FCS | 14 |
| Siemens | | Sirona Operatory | | 12V,50W | BW.64440 | 12 |
| Siemens | | X-Ray | | 12V,20W | BW.64115 | 55 |
| Simplified Systems | | Communicator | | 28V,.04A | BW.85 | 40 |
| Sirona | | LED High Speed Fiber Optic Handpiece | | 6 Pin Coupler | BW.VSBL | - |
| Solarc | | Patterson Dental | AL1824, WA09500, 09500 | 24W Solarc | BW.09500-USH | 8 |
| Solid State Systems | | Initiator | 6220002-290 | 15V,150W | BW.EFR | 1 |
| Sony | | 20/20 Carry Cam | | 21V,150W | BW.EJA | 1 |
| S.S. White | | Profile VCS | 17244 | 12V,100W | BW.JCRM12V-100W | 4 |
| Spectraview | | Intra & Extraoral Camera | | 12V,75W | BW.EFN | 1 |
| Springhealth Penn-Med | | Original Wand | | 120V,150W | BW.ESD | 1 |
| Springhealth Penn-Med | | Powerlite 75; 100; 310 Gun | BU75 | 12V,75W | BW.JCRM12V-75W | 4 |
| Springhealth Penn-Med | | 200L Wand | | 21V,150W | BW.EJV | 1 |
| Star Dental | | Handpiece for Quick Connect Swivel | 260986 | 4.5-5.5V 2 Piece Pack | BW.60986 | 32 |
| Star Dental | | Starflite | 61351 | 21V,150W | BW.DNF | 6 |
| Star X-Ray | | DE 94; 98; 100; 120; 150; 160; 200; SW100 Viewbox | | 14W Fluorescent | BW.F14T8.D | 62 |
| Star X-Ray | | DE 99 | | 15W Fluorescent | BW.F15T8.D | 62 |
| Star X-Ray | | Endodontic Viewer | | 8W Fluorescent | BW.F8T5.D | 61 |
| Star X-Ray | | Mini Viewer | | 6W Fluorescent | BW.F6T5.D | 61 |
| Star X-Ray | | Starcure 350 EH | | 120V,150W | BW.JCR120V150WB | 1 |
| Star X-Ray | | Starcure 360 | | 9W Blue | BW.CF9DS.BLUE | 59 |
| Star X-Ray | TR 130 Cephalometric | | 8W Fluorescent | BW.F8T5.D | 61 | |
| Star X-Ray | 612 Duplicator | | 120V,7W Orange | BW.7C7-ORANGE | 25 | |
| Steris | Examiner 10 | | 120V,150W | BW.ESP | 19 | |
| Summit Dental Sys. | Overhead Light | | 24V,70W | BW.JA24V-70WD | 28 | |

▲ Multiple Brands Other Than Brand Shown May Be Available – See Price List

Call Bulbworks • 800-334-2852 (973-584-7171) For Any Manufacturer Or Model Not Listed!

E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com EUPS2020A DENT

Equipment To Lamp Cross Reference Guide
For Pricing See Separate Bulbworks Price Sheet

| | Manufacturer | Model | Mfr's Part No. | Description | Bulbworks Part No.▲ | Fig. |
|----------|----------------------|---|---------------------------|--------------------------------|---|----------------|
| T | Teledyne Getz | Focus Activator | 4501 | 12V,100W | BW.JCRM12V-100W | 4 |
| | Theta Communications | Five-Lite – Square Lenses 1972 and before | | 6.3V,.15A | BW.40 | - |
| | Theta Communications | Five-Lite – Square Lenses 1973 and after | | 6.3V,.15A | BW.47 | 54 |
| | Theta Communications | Five-Lite – Square Lenses 1997 and before | | 6.3V,.20A | BW.381 | 53 |
| | Theta Communications | Five-Lite – Square Lenses 1998 and after | | 6.3V,.15A | BW.47 | 54 |
| | Theta Communications | Ten-Lite – Square Lenses 1997 and before | | 28V,.06A | BW.376 | 53 |
| | Theta Communications | Ten-Lite – Square Lenses 1998 and after | | 28V,.04A | BW.1819 | 54 |
| | Theta Communications | Ten-Lite – Round Lenses | | 28V,.04A | BW.1819 | 54 |
| | Theta Communications | Status System – Round Lenses of Units with 600 Part Numbers | | 6.3V,.15A | BW.47 | 54 |
| | Theta Communications | Status System – Round Lenses of Units with 700 Part Numbers | | 28V,.04A | BW.1819 | 54 |
| | 3M | Visilux I; II | 78-8045-4571-9 | 12V,75W | BW.JCRM12V-75W | 4 |
| | 3M | Visilux 5510 | | 12V,50W | BW.JC12V-50W | 12 |
| | 3M | XL1500 | 3367-5267 | 10V,52W | BW.3367-5267* | 3 |
| | 3M | XL2500; 3000 | 78-8045-4571-9 | 12V,75W | BW.JCRM12V-75W | 4 |
| | 3M-ESPE | Elipar I; Elipar Visio Beta | 70-2011-0200-4 | 13V,75W | BW.70201102004 | 2 |
| | 3M-ESPE | Elipar II; Elipar Mini Spot High Light | 70-2011-0199-8 | 12V,75W | BW.64613-IMP | 4 |
| | 3M-ESPE | 2500 | 78-8131-1662-7 | 12V,75W | BW.64613-IMP | 4 |
| | Trojan Research | Ultraeye | | 21V,150W | BW.EKE | 1 |
| | Trophy | Intra Oral | | 12V,75W | BW.EFN | 1 |
| | U | Unitek/3M | Ortholux; XT | 78-8045-4571-9 | 12V,75W | BW.JCRM12V-75W |
| V | Vacudent | | | 24V,150W | BW.FDS.DZE | 17 |
| | Vector | VKS6 Handpiece | VRB Bulb | Generic 1 Pack | BW.VRB-GEN1 | 41 |
| | Vector | VKS6 Handpiece | VRB Bulb | Generic 3 Pack | BW.VRB-GEN3 | 41 |
| | Vidicam | | | 120V,150W | BW.JCR120V150WB | 1 |
| | Vivadent | Astralis V, VII | | 12V,75W | BW.14552 | 3 |
| | Vivadent | Astralis 10 | 561398 | 13V,100W | BW.561398 | 3 |
| | Vivadent | Heliolux; Heliolux II | 7412, 133762 | 14V,35W | BW.13165-IMP | 3 |
| | Vivadent | Heliolux DLX; GT | 7512 | 12V,75W | BW.JCRM12V-75W | 4 |
| | Vivadent | Heliomat B Single | 7210 | 21V,150W | BW.EJA | 1 |
| | Vivadent | Heliomat H10 Single; Heliomat Multi | 7110 | 15V,150W | BW.EFR | 1 |
| W | W & H | W & H Handpiece (Can Use Standard Bulb or LED) | 22076 | LED Version - Rated as 3.5V | BW.22076-LED Can Be Used To Replace Standard 22076 Bulb | 45 |
| | Weber | Panavision | | 120V,125W | BW.125T10.P | 20 |
| | Weber | SK81G | | 12V,100W | BW.JA12V-100WD | 28 |
| | Weber | SK100 Skylight | | 12V,100W | BW.FDT | 18 |
| | Weber | SK800; SK900G Skylight | SK580 | 14.5V,90W | BW.EPV | 1 |
| | Weil Dental | Elcelux | | 14V,35W | BW.13165-IMP | 3 |
| | Welch Allyn | DA-L-400 | WA09800 | 21W | BW.WA09800 | 8 |
| | Welch Allyn | Denlite Chair Site Mirror | DP1020 | 1 Piece Pack | BW.DP1020 | 36 |
| | Welch Allyn | Denlite Portable Mirror | DP1520 | 1 Piece Pack | BW.DP1520 | 36 |
| | Welch Allyn | Exam Light II | WA04200 | 12V,35W | BW.EPN | 1 |
| | Welch Allyn | Solarc, Reveal DLS400 | AL1824, WA09500, 09500 | 24W Solarc | BW.09500-USH | 8 |
| | S.S.White | Profile VCS | 17244 | 12V,100W | BW.JCRM12V-100W | 4 |
| | Wykle Research | Witelite 1500 | 3052 | 12V,75W | BW.JCRM12V-75W | 4 |

▲ Multiple Brands Other Than Brand Shown May Be Available – See Price List

* Same as Caulk/Dentsply #644751 operated at different voltage.

Call Bulbworks • 800-334-2852 (973-584-7171) For Any Manufacturer Or Model Not Listed!
E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com

**Replacement Lamps For Dental Equipment
Distributor/Industry Number Cross Reference Guide**

| Distrib. No. (Efos) | Distrib. No. (Efos & Caulk) | Distrib. No. (Chapman-Huffman) | Distrib. No. (Quala) | Philips No. | Osram No. | GE, Industry Std., Or ANSI No. | Bulbworks No. |
|------------------------|--------------------------------|-----------------------------------|-------------------------|-------------|-----------|-----------------------------------|-----------------|
| 3008 | 640 30080 | 16-432-00 | - | - | - | EZK | BW.JCR120V150WB |
| 3009 | 644751 | 16-433-00 | QUA-230 | 13298 | - | - | BW.3367-5267 |
| 3010 | 640 30100 | - | - | - | 64613 | 64613 | BW.64613 |
| 3011 | 640 30110 | 16-360-00 | QUA-160 | 5995 | - | EJM | BW.EJM |
| 3014 | 640 30140 | 16-330-00 | QUA-165 | 6423 | 64634 | EFR | BW.EFR |
| 3015 | 640 30150 | 16-410-00 | QUA-170 | - | - | ETJ | BW.ETJ |
| 3016 | 640 30160 | 16-420-00 | QUA-175 | 13865 | 64617 | JCRM12V-75W | BW.JCRM12V-75W |
| 3017 | 640 30170 | 16-350-00 | QUA-180 | 13164 | 64644 | EJL | BW.EJL |
| 3018 | 640 30180 | 16-340-00 | QUA-185 | - | - | EJA | BW.EJA |
| 3019 | 640 30190 | 16-300-00 | - | - | - | DNE | BW.DNE |
| 3020 | 640 30200 | 16-320-00 | QUA-190 | 6834 | 64627 | EFP | BW.EFP |
| 3021 | 640 30210 | 16-430-00 | QUA-195 | - | 64624 | JCRM12V-100W | BW.JCRM12V-100W |
| 3022 | 640 30220 | 16-370-00 | QUA-200 | - | - | EJV | BW.EJV |
| 3023 | 640 30230 | - | - | - | - | DNF | BW.DNF |
| 3024 | 640 30240 | 16-390-00 | - | - | - | ELZ | BW.ELZ |
| 3025 | 640 30250 | - | - | - | - | - | BW.650197 |
| 3026 | 644725 | 16-431-00 | QUA-205 | 13165 | - | JCRM14V-35W | BW.13165 |
| 3027 | 640 30270 | 16-380-00 | QUA-210 | - | - | EKE | BW.EKE |
| 3028 | 640 30280 | 16-400-00 | QUA-215 | - | - | ESD | BW.ESD |
| 3031 | 640 30310 | 16-054-00 | - | - | - | 041-317-00 | BW.041-317-00 |
| 3032 | 640 30320 | - | - | - | 64440 | JC12V-50W | BW.JC12V-50W |
| 3033 | 640 30330 | - | - | - | - | - | *** |
| 3034 | 640 30340 | 16-445-00 | - | - | - | ESY | BW.ESY |
| 3035 | 640 30350 | 16-436-00 | QUA-220 | 6853 | 64615 | EFN | BW.EFN |
| 3036 | - | - | - | - | - | Q50MR16/FL/1 - ENL | BW.ENL |
| 3038 | 640 30380 | 16-437-00 | - | 7748 | 64655 | EHJ | BW.EHJ |
| 3040 | 640 30400 | 16-435-00 | - | 6847 | 64607 | EFM | BW.EFM |
| 3047 | 640 30470 | - | - | - | - | FC8T9/CW | BW.FC8T9.CW |
| 3048 | 640 30480 | - | - | - | - | F6T5/D | BW.F6T5.D |
| 3049 | 640 30490 | - | QUA-140 | - | - | F8T5/D | BW.F8T5.D |
| 3049B | 640 3049B | - | - | - | - | F8T5.BLB | BW.F8T5.BLB |
| 3050 | 640 30500 | - | - | - | - | F6T5/CW | BW.F6T5.CW |
| 3051 | 640 30510 | - | - | - | - | F8T5/CW | BW.F8T5.CW |
| 3052 | 640 30520 | - | - | - | - | F13T5/CW | BW.F13T5.CW |
| 3053 | 640 30530 | - | QUA-155 | - | - | F14T12/CW | BW.F14T12.CW |

*** Not Available Through Bulbworks As Of April, 2004

Continued On The Next Page...

**For Additional Cross Reference Information Contact
Bulbworks, Inc. • P.O. Box 586, Succasunna, NJ 07876
CALL: 800-334-2852 (973-584-7171) • FAX: 800-335-2858 (973-584-0300)
E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com**

Replacement Lamps For Dental Equipment Distributor/Industry Number Cross Reference Guide

Continued From The Previous Page

| Distrib. No. (Efos) | Distrib. No. (Efos & Caulk) | Distrib. No. (Chapman-Huffman) | Distrib. No. (Quala) | Philips No. | Osram No. | GE, Industry Std., Or ANSI No. | Bulbworks No. |
|------------------------|--------------------------------|-----------------------------------|-------------------------|-------------|------------|-----------------------------------|-----------------|
| 3054 | 640 30540 | - | QUA-150 | - | - | F14T12/D | BW.F14T12.D |
| 3055 | 640 30550 | - | - | - | - | F14T8/CW | BW.F14T8.CW |
| 3056 | 640 30560 | - | QUA-145 | - | - | F14T8/D | BW.F14T8.D |
| 3057 | 640 30570 | - | - | - | - | F15T8/CW | BW.F15T8.CW |
| 3058 | 640 30580 | - | - | - | - | F15T8/D | BW.F15T8.D |
| 3059 | 640 30590 | 16-411-00 | - | 7724 | 64623 | EVA | BW.EVA |
| 4070 | 640 40700 | 16-020-00 | QUA-115 | 7158 | 64640 | JC24V-150W | BW.FCS |
| 4071 | 640 40710 | 16-040-00 | QUA-100 | 12336 | 64151 | JA12V-55WD (H3) | BW.JA12V-55WD |
| 4072 | 640 40720 | 16-035-00 | QUA-105 | 125T10/P | - | 125T10/P | BW.125T10.P |
| 4073 | 640 40730 | 16-015-00 | - | - | - | Q250/MC/CL-120V EHT | BW.EHT |
| 4074 | 640 40740 | 16-025-00 | QUA-110 | 5974 | 64643 | FDS/DZE | BW.FDS.DZE |
| 4075 | 640 40700 | 16-020-00 | QUA-115 | 7158 | 64640 | JC24V-150W | BW.FCS |
| 4076 | 640 40760 | 16-050-00 | QUA-120 | - | - | Q150T4/CL-25V | BW.Q150T4.CL25V |
| 4077 | 640 40770 | 16-050-00 | QUA-120 | - | - | Q150T4/CL-25V - GE Brand | BW.Q150T4CL25GE |
| 4078 | 640 40780 | - | - | - | - | - | BW.PC TOOL |
| 4079 | 640 40790 | 16-045-00 | QUA-125 | - | - | Q75CL/DC-28V | BW.Q75CL.DC-28V |
| 4080 | 640 40800 | 16-010-00 | - | - | - | EPV | BW.EPV |
| 4081 | 640 40810 | 16-030-00 | - | 5973 | 64628 | FDT | BW.FDT |
| 4082 | - | - | - | - | - | 100A21/SP | BW.100A21SP |
| 4084 | 640 40840 | 16-351-00 | - | - | - | ENZ | BW.ENZ |
| 4085 | 640 40850 | 16-011-00 | - | 6642 | - | ESX | BW.ESX |
| 5616 | - | 16-320-00 | QUA-190 | 6834 | 64627 | EFP | BW.EFP |
| 8023 | - | - | - | - | - | FC8T9/CW | BW.FC8T9.CW |
| 96000 | 640 96000 | - | - | - | CF9DS/BLUE | CF9DS/BLUE | BW.CF9DS.BLUE |
| - | - | 16-012-00 | - | - | - | JCD24V-75W | BW.JCD24V-75W |
| - | - | 16-058-00 | - | - | - | 650064 | BW.650064-GEN2 |
| - | - | 16-061-00 | QUA-135 | - | - | 553-3881 | BW.553-3881 |
| - | - | 16-062-00 | - | - | - | Y900529 | BW.Y900529-GEN1 |
| - | - | 16-063-00 | - | - | - | 42001021 | BW.42001021GEN1 |
| - | - | 16-440-00 | QUA-225 | - | - | 21237 | BW.21237 |
| - | - | - | QUA-130 | - | 64451 | - | BW.64451 |
| - | - | - | QUA-131 | 14623 | - | - | BW.041-709-00 |

For Additional Cross Reference Information, Contact:

Bulbworks, Inc. • P.O. Box 586, Succasunna, NJ 07876
CALL: 800-334-2852 (973-584-7171) • FAX: 800-335-2858 (973-584-0300)
E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com

Getting To Know Us... (Terms & Conditions)

ORDERING IS EASY...

Call Us TOLL FREE • 800-334-BULB(2852) • Our business hours are 8:30 AM to 5:00 PM Eastern Time Monday through Thursday and 8:30 AM to 4:00 PM on Friday. Outside of business hours? Use our after hours message machine and your call will be returned during the next business day or you may...

Fax Us ON OUR 24 HOUR TOLL FREE FAX LINE • 800-335-2858 • All orders received by fax are acknowledged promptly either by phone or return fax.

Electronic Orders • Send orders to orders@bulbworks.com • All orders placed by e-mail will be acknowledged promptly by return e-mail.

Order Online • www.bulbworks.com • All orders placed online will be acknowledged promptly by return e-mail.

No Minimum Order • We appreciate your business no matter how small or large your order may be. Orders with a merchandise value of under \$10.00 will be subject to a Small Order Handling Charge of \$3.00 (In addition to Merchandise and Shipping).

Fast Shipment • Orders received by 3:00 PM Eastern Time will be **shipped the same day** that we receive them. Orders are normally shipped by UPS but all express services are available. Bulbworks carefully inspects all lamps prior to shipment and we use the best packaging methods available.

Terms • We accept MasterCard, Visa, Discover, & American Express. Our terms for *approved* open accounts are 2% 10 Days, Net 30, freight prepaid and added to the invoice. We would be glad to open an account for your convenience. Just ask!



Order With Confidence • All items carry a full warranty. Bulbworks sells only top quality brand name products - the same brands used by the equipment manufacturers. Should product failure, breakage, or defect occur, contact us for immediate replacement!

Original Manufacturer's Packaging • Each lamp is shipped in its original packaging from the lamp manufacturer – Eiko, Philips, GE, Ushio, Osram - Sylvania, Welch Allyn, and others. Part numbers are clearly labeled and are easy to read to eliminate confusion at reorder time. Helpful equipment reference information appears on every carton.

All Prices Are Subject To Change Without Notice. Some Prices May Be For Stock On Hand Only.
No Part Of This Catalog May Be Reproduced Without Permission.



P.O. Box 586 Succasunna, NJ 07876
CALL: 800-334-BULB (2852) • (973-584-7171)
FAX: 800-335-2858 • (973-584-0300)
E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com



**Need Something Besides Dental Bulbs? Stop Looking & Call Us!
We Have What You Need! We Stock All Types Of
Bulbs, Sockets, and Batteries Including:**

- **MICROSCOPE & SCIENTIFIC BULBS** Replacement Lamps For American Optical, Bausch & Lomb, Jena, Leitz, Nikon, Olympus, Swift, Unitron, Zeiss, And Others. Lamps For Biological, Biomedical, Educational, Forensic, Industrial, Metallurgical, Research, And Other Applications. Fluorescent And LED Microscope Illuminators.
- **MEDICAL BULBS** Replacement Lamps For Operatory, Surgical, Clinical, Fiber Optic, Laser, Analytical, Veterinary, X-Ray Viewing, And Other Applications. Lamps For All Types Of Medical Equipment.
- **INSPECTION** Replacement Lamps For Optical Comparators Including Baty Gagemaster, J & L, Micro Vu, Mitutoyo, Nikon, OGP, Scherr Tumico, Vermont Precision, And Others. Replacement Lamps for Microscopes And Other Inspection Devices.
- **GRAPHIC ARTS BULBS** Replacement Lamps For AB Dick, Agfa, Amergraph, Brown, Burgess, Douthitt, Itek, Kenro, NuArc, Olec, Sandmar, Screen Theimer, 3M, And More.
- **OPTICAL AND OPHTHALMIC BULBS** Replacement Lamps For American Optical, Bausch & Lomb, Copeland, Haag Streit, Marco, Mentor, Nikon, Reichert, Topcon, Welch Allyn, And Others. Lamps For Optical And Ophthalmic Equipment Including Lensometers, Ophthalmoscopes, Projectors, Retinoscopes, Slit Lamps, Tonometers, And Others.
- **GENERAL AND SPECIALTY LIGHTING** Replacement Lamps And LEDs For Color Enhancement And Correction, Commercial, Decorative, Energy Saving, General Service, Long Life, Office, Problem Solving, Rough Service, And Other Applications.
- **GERMICIDAL, BLACK LIGHT, AND OTHER UV BULBS** Replacement Lamps For Air And Water Purification. For Use In Hospitals, Laboratories, Manufacturing, Food Handling, And Many Other Applications.
- **PROJECTOR AND AUDIO VISUAL BULBS** Replacement Lamps For Overhead Projectors, Slide And Movie Projectors, Microfilm and Microfiche Readers, Micrographic And Other Applications.
- **SOCKETS AND LAMP HOLDERS** Replacement Sockets For All Lamp Types Including Halogen, Incandescent, Vacuum, & Others. Sockets From Bender & Wirth, Buhl, Gilway, Christiana Industries, Cole Hersee, Gilway, Sylvania, Ushio, And Others.
- **BATTERIES** Batteries For All Types Of Devices. Alkaline, Lithium, Nickel Metal Hydride, Heavy Duty, Zinc Air, Rechargeable, Industrial, Hearing Aid, Photo, Electronic, Medical, Keyless Entry, Watch and Others. All Sizes Available.

We Supply All Specialty Lamp Types Including: Fluorescent, Incandescent, Halogen, Mercury Vapor, Metal Halide, Miniature, Quartz, Xenon, And Other Specialty Bulbs.

BULBWORKS Supplies Only Top Quality Brand Names - The Same Brands Used By The Equipment Manufacturers. Lamps From Eiko, GE, Hosobuchi, Narva, Nikon, Opti Quip, Osram Sylvania, Philips, Ushio, Welch Allyn, ORC, Voltarc, Zeiss, And Many Other Specialty Lamp Manufacturers

Call BULBWORKS For Your Next Replacement Lamp Requirement!

Bulbworks, Inc. • P.O. Box 586 • Succasunna, NJ 07876

CALL: 800-334-BULB(2852) • (973-584-7171)

FAX: 800-335-2858 • (973-584-0300)

E-mail: bulbwork@bulbworks.com • Order Online: www.bulbworks.com